INSTRUCTIONS FOR USE FOR:

en

English

dk

Dansk

nl

Nederlands

fi

Suomi

fr

Français

de

Deutsch

gr

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it

Italiano

no

Norsk

pt

Português

es

Español

se


Svensk

INSTRUCTIONS FOR USE: GORE INTRODUCER SHEATH WITH HEMOSTASIS VALVE

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions. Failure to do so may result in complications.

DESCRIPTION

The GORE Introducer Sheath with Hemostasis Valve consists of an outer sheath with hemostasis valve and sideport, as well as an inner dilator with a tapered tip. The main outer sheath tubing is connected at the proximal (trailing) end to a hemostasis valve with side port tubing that is connected to a plastic 3-way stopcock valve. The side port is used for flushing the sheath. The sheath is introduced into the vascular system with the aid of the inner dilator which has a luer lock at the proximal (trailing) end and a tapered tip at the distal (leading) end. The hemostasis valve at the proximal (trailing) end of the sheath conforms and seals around guidewires and catheters to reduce blood leakage from the sheath.

	Catalogue Number		
	PG123000	PG183000	PG203000
A	37 cm	40 cm	40 cm
B	2.0 cm	5.1 cm	5.1 cm
C	12 Fr (4 mm)	18 Fr (6 mm)	20 Fr (7 mm)
D	.040"	.037"	.037"
E	4.7 mm	6.8 mm	7.6 mm
F	30 cm	30 cm	30 cm

NOTE: Dimensions listed are nominal

PACKAGE CONTENTS

- Sheath
- Dilator

INTENDED USE

The GORE Introducer Sheath is used to facilitate the introduction of guidewires, catheters and other interventional medical devices into the vascular system, and to minimize blood loss associated with such introduction.

PRECAUTIONS

- Examine packaging and device before use. Do not use if either the packaging or device is damaged or if sterile barrier has been compromised.
- The GORE Introducer Sheath is supplied sterile. The package will serve as an effective barrier until the "use by" (expiration) date printed on the box. Do not resterilize; for single use only.
- Store in a cool dry place.
- Individual patient anatomy and physician technique may require procedural variations.
- Verify sheath, device, catheter and accessory components size compatibility prior to use.

WARNINGS

- Do not alter this device. Alterations may impair function.
- Do not attempt to use a guidewire with a maximum diameter greater than specified on package label.
- This device should only be used by physicians thoroughly trained in the technique of catheter delivery systems.
- The GORE Introducer Sheath should only be advanced or retracted under fluoroscopic guidance, and should only be advanced with the dilator fully inserted.
- Do not attempt to advance or withdraw guidewire, catheter, or other interventional medical device through introducer sheath and/or dilator if resistance is felt. Use fluoroscopy to determine cause. Continued advancement or retraction against resistance may result in damage to the vessel, or breakage of the guidewire, catheter or interventional medical device. Continued advancement or retraction against resistance may result in serious injury.
- Advance dilator/sheath assembly with a twisting motion to avoid damage to the sheath or vessel.
- Do not attempt to insert a catheter or other interventional medical device having a diameter larger than the introducer sheath size indicated. Device damage or breakage may occur.

DIRECTIONS FOR USE

Sheath Preparation

1. Verify proper size (diameter and length) is selected.
2. Remove the Introducer Sheath from its packaging and examine for possible damage or defects. Do not use any damaged or defective devices.
3. Flush the dilator, introducer sheath, and sideport with heparinized intravenous fluid.
4. When the sheath will remain in a vessel for an extended period, consider using a continuous drip of heparinized intravenous fluid under pressure administered through the sideport connection.

Sheath Use

1. Insert the dilator tip through the valve and completely into the sheath until the dilator hub comes in contact with the hemostasis valve. This ensures that the tapered portion of the dilator is beyond the end of the sheath.
NOTE: When advancing the sheath, ensure that the dilator remains fully inserted into the sheath and that the dilator hub remains in contact with the hemostasis valve.
2. Follow normal accepted practice for vessel puncture or incision and guidewire insertion.
3. Using fluoroscopic guidance, advance the dilator/sheath over the guidewire as a unit; do not allow dilator to back out of the (separate) sheath while advancing. Stop advancement of the assembly if there is resistance. Investigate the cause of resistance before proceeding. Carefully advance the assembly until it is at the desired location. Advance dilator/sheath assembly with a twisting motion to avoid damage to the sheath or vessel.
4. Hold the sheath steady and maintain the guidewire position while withdrawing the dilator from the sheath, over the guidewire until it is completely removed from the guidewire.
5. While maintaining the position of the guidewire, advance the selected catheter or other interventional medical device over the guidewire into the sheath, taking care to keep the sheath assembly as straight as possible outside the body and avoid kinking.
6. When exchanging different catheters and devices through the introducer sheath care should be taken to maintain proper guidewire and sheath positions within the vascular system.
7. Carefully support all wires, catheters and devices while pushing across the hemostasis valve.
8. Upon removal of the sheath, precautions should be taken to prevent bleeding, vessel damage, or other serious injury.

CAUTIONS

- Advance and remove sheath only under fluoroscopic guidance.
- Do not attempt sheath advancement without guidewire and dilator secured in place. Severe vascular damage and/or injury may occur.
- Insertion into and removal from artery may cause excessive bleeding and/or other complications.
- Do not continue advancement or retraction of the catheter or other interventional medical device into and out of the introducer if there is resistance. Determine the cause of resistance before proceeding.

DEFINITIONS

 Use By

 Attention, See Instructions for Use

 Do Not Reuse

REF Catalogue Number

LOT Batch Code

EU REP European Authorized Representative

STERILE

Contents sterile unless package has been opened or damaged.

STERILE EO

Contents sterile unless enclosed package has been opened or damaged. Sterilized by ethylene oxide.

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© 2002, 2003 W. L. Gore & Associates, Inc.
GORE and designs are trademarks of W. L. Gore & Associates.
W. L. Gore & Associates, Inc.
W. L. Gore & Associés, S.A.R.L.