


PROVEN PATENCY.
MEASURABLE VALUE.




Strength of leadership



GORE® PROPATEN® Vascular Graft


A leading prosthetic vascular graft for lower extremity revascularization, specifically designed to reduce the risk of acute graft thrombotic failure. With more than a decade of strong performance that includes improving outcomes and reducing interventions, this longstanding bypass graft helps deliver both proven clinical and economic value for patients and hospitals.

See the proof at goremedical.com/eu/propaten



600K+

Devices
implanted



30+

Clinical studies
over 10 years



2,000+

Limbs
studied




17+

Years of
performance

A decade of performance


GORE® PROPATEN® Vascular Graft

Above-knee bypass primary patency.



GORE® PROPATEN® Vascular Graft

Below-knee bypass primary patency.



* Overall weighted average primary patency is based on data from 15 peer-reviewed publications meeting pre-determined inclusion criteria. Visit propatenperformance.com to see inclusion criteria, explore the data, and see publications.

† Data not reported.

Study size (N) reflects the initial cohort size of the study.

Proven patency

By substantially reducing acute graft thrombosis within hours after implantation, the CBAS® Heparin Surface on the GORE® PROPATEN® Vascular Graft provides clinical benefits that standard ePTFE grafts do not.¹



Fewer occlusions

50% reduction in risk of graft occlusion compared to standard ePTFE in critical limb ischemia (CLI) patients.²




Improved patient outcomes

Higher primary and secondary patency and higher limb salvage for below-knee bypass compared to standard ePTFE from 1–3 years.³


Improved clinical outcomes

Primary patency




GORE® PROPATEN® Vascular Graft³ Standard ePTFE³

Secondary patency




Limb salvage




Below-knee cumulative amputation-free survival


GORE® PROPATEN® Vascular Graft Standard ePTFE



GORE® PROPATEN®
Vascular Graft
long-term value*



Improvement in
amputation-free survival



Decrease in
revision procedures

Amputation-free survival (avoided loss of limb or life) is the average reported mortality rate for standard ePTFE and the average reported amputation rates for standard ePTFE and GORE® PROPATEN® Vascular Graft.⁴

* Based on the three-year published clinical performance and economic model.

Lasting thromboresistance. Proven technology.*

CBAS® Heparin Surface

The CBAS® Heparin Surface of the GORE® PROPATEN® Vascular Graft consists of a proprietary covalent end-point bond that preserves the active site, thus retaining heparin's anticoagulant activity.



Proven heparin availability

Performance-ready heparin active site.^{4,6}



Proven heparin bioactivity

Unmatched, persistent ability to take up antithrombin.^{1,5}



Proven lasting thromboresistance

Improved surface hemocompatibility resulting from heparin availability and bioactivity.^{1,4-7}


In vivo canine carotid artery interposition model⁶

The CBAS® Heparin Surface of a 3 mm diameter GORE® PROPATEN® Vascular Graft (top) remains free of thrombus, while the 3 mm diameter control ePTFE graft (bottom) is covered with thrombus in an acute two-hour in vivo canine carotid artery interposition model.⁶

GORE® PROPATEN® Vascular Graft




Standard ePTFE



* CBAS® Heparin Surface. W. L. Gore & Associates Web site. <https://www.goremedical.com/eu/cbas/references>. Accessed September 25, 2019.

Sustained heparin bioactivity^{5,8}




Eight years (Explant after 2,939 days)

Heparin bioactivity detected above the level required for thromboresistance in an eight-year human explant. No adherent thrombus was found.⁸

- Femoral to posterior tibial bypass with polyester Linton patch
- Distal anastomosis occluded

Proprietary end-point covalent bonding




The anticoagulant function of heparin is dependent on the bioavailability of an active site within the molecule. Some methods of covalent heparin bonding damage and/or obstruct the active site and hence destroy heparin's anticoagulant activity.

Visit goremedical.com/eu/cbas to learn more

Features Gore's CBAS® Heparin Surface, the proven heparin bonding technology for lasting thromboresistance, used in many of Gore's interventional and vascular surgery products. End-point covalent bonding keeps heparin anchored to the graft surface, while the bioactive site remains free to interact with the blood to help prevent clotting.⁷

References

1. Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® BioActive Surface heparin immobilization. *European Journal of Vascular & Endovascular Surgery* 2003;25(5):432-437.
2. Lindholt JS, Gottschalksen B, Johannesen N, et al. The Scandinavian PROPATEN® Trial – 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses – a randomized clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673.
3. GORE® PROPATEN® Vascular Graft. W. L. Gore & Associates Web site. <https://www.goremedical.com/products/propaten/references>. Accessed May 15, 2019.
4. Gore S, Andersson J, Biran R, Underwood C, Riesenfeld J. Heparin surfaces: impact of immobilization chemistry on hemocompatibility and protein adsorption. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 2014;102(8):1817-1824.
5. Freeman J, Chen A, Weinberg RJ, Okada T, Chen C, Lin PH. Sustained thromboresistant bioactivity with reduced intimal hyperplasia of heparin-bonded PTFE Propaten Graft in a chronic canine femoral artery bypass model. *Annals of Vascular Surgery* 2018;49:295-303.
6. Biran R, Pond D. Heparin coatings for improving blood compatibility of medical devices. *Advanced Drug Delivery Reviews* 2017;112:12-23.
7. CBAS Heparin Surface. W. L. Gore & Associates Web site. <https://www.goremedical.com/eu/cbas/references>. Accessed September 25, 2019.
8. Biran R. *Heparin Activity and Concentration Values from 52 and 96 Month Duration Clinical Explants of the GORE® PROPATEN® Vascular Graft*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2011. [Work plan]. WP104647.



 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

BAYER and ASPIRIN are trademarks of Bayer Intellectual Property GmbH or its affiliated companies.

MEDLINE is a trademark of Medline Industries, Inc.

SANOFI-AVENTIS and PLAVIX are trademarks of Sanofi-Aventis.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

GORE, GORE-TEX, *Together, improving life*, PROPATEN and designs are trademarks of W. L. Gore & Associates.

© 2022 W. L. Gore & Associates GmbH 22758838-EN OCTOBER 2022

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673

United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

