


Gore VIPER Clinical Study One-Year Results



© 2014 W. L. Gore & Associates, Inc.

Product with radiopaque markers planned for European availability in 2016.

GORE®, VIABAHN®, and designs are trademarks of W. L. Gore & Associates. © 2007-2015 W. L. Gore & Associates, Inc.



Gore VIPER Clinical Study Centers and Principal Investigators


North County Medical Center, Oceanside, California	Richard Saxon, MD Study Principal Investigator
Boone Hospital Center, Columbia, Missouri	John Adams, MD
Riverside Methodist Hospital, Columbus, Ohio	Gary Ansel, MD
St. Luke's Medical Center, Milwaukee, Wisconsin	Tanvir Bajwa, MD
Heritage Valley Health System, Beaver, Pennsylvania	Richard Begg, MD
Vascular Surgical Associates, PC, Austell, Georgia	Arun Chervu, MD
University of California, San Francisco, California	Charles Eichler, MD
Baylor University, Dallas, Texas	Dennis Gable, MD
Mercy Hospital, Chicago, Illinois	Paul Jones, MD
Holy Cross Hospital, Fort Lauderdale, Florida	Michael Rush, MD
St. Elizabeth's Medical Center, Boston, Massachusetts	Peter Soukas, MD

Gore VIPER Clinical Study Overview

GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface for treatment of long SFA disease

Objective	Evaluate the performance of GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface (W. L. Gore & Associates, Inc.) in treating long-segment SFA disease (> 5 cm in length)
Design	Single-arm, prospective, 12 sites, 120 limbs
Primary Endpoints	<p>Primary patency at 12 months</p> <ul style="list-style-type: none">• No evidence of restenosis or occlusion within the originally treated lesion based on CDUS; PSVR < 2.5;• No angiographic evidence of stenosis > 50% if CDUS is uninterpretable or unavailable <p>Proportion of subjects experiencing major procedure related adverse events within 30 days of procedure</p>
Secondary Endpoints	<p>Primary assisted patency</p> <p>Secondary patency</p> <p>Device-related major adverse events at 12 months</p>

Endoprosthesis Description



Product with radiopaque markers planned for European availability in 2016.

CBAS[®] is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.


GORE[®], VIABAHN[®], and designs are trademarks of W. L. Gore & Associates. © 2007-2015 W. L. Gore & Associates, Inc.

Contoured Proximal Edge

Manufacturing change during study

IVUS demonstrates device apposition to artery in canine model

Post-mortem dissection demonstrates device apposition to artery




Non-Contoured



Contoured

CBAS[®] Heparin Surface

- End-point covalent bonding
 - CBAS[®] Heparin Surface technology allows for retention of bioactivity.
- Sustained bioactivity¹
 - Active site catalytically facilitates antithrombin-thrombin complex formation and then becomes available to repeat the reaction.
- Intended to provide a thromboresistant surface
 - Clinical history: Long-term activity and safety.²



1. P.C. Begovac, R.C. Thomson, J.L. Fisher, A. Hughson, A. Gallhagen. Improvements in GORE-TEX[®] vascular graft performance by Carmeda[®] bioactive surface heparin immobilization. *European Journal of Vascular & Endovascular Surgery* 2003; 25(5): 432-437.
2. Lindholt JS, Gottschalksen B, Johannesen N, *et al.* The Scandinavian Propaten[®] Trial – 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses – a randomised clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673.


Lesion Characteristics

	Gore VIPER Clinical Study
Limbs Enrolled	119
Treated Occlusions	56%
Lesion Length	19 cm
Lesion Calcification	
none-mild	39%
moderate-severe	61%
Tibial Runoff	
One vessel	21%
Two vessel	33%
Three vessel	46%
TASC II Lesion Classification	
Type A	14%
Type B	25%
Type C	29%
Type D	31%

Safety – Major Adverse Events

- **Primary Endpoint: 30-day procedure-related MAE**
 - One event, (0.8%): surgical bypass after target lesion occlusion
- **Secondary Endpoint: One-year device-related MAE**
 - Zero events
- **Major Adverse Events (MAE): require significant therapy, including unplanned increase in the level of care, permanent sequelae, hospitalization, or death. Repeat interventions, stenosis, and occlusions are not adverse events.**

One-Year Patency




103 / 119 limbs available for follow-up at 12 months

One-Year Primary Patency by Subgroup


	Primary Patency
Overall	73%
Device Diameter	
5 mm (n= 23)	79%
6 mm (n= 85)	69%
7 mm (n= 8)	100%
Lesion Length	
≤ 20 cm (n= 68)	75%
> 20 cm (n= 51)	70%

Effects of Device Sizing: Proximal



Device oversizing assessed by independent Core Lab, data on file

Oversizing




Conclusions

- **The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface Exhibits 73% Patency in Long SFA Lesions**
 - Patency is independent of lesion length
 - Long lesions (> 20 cm) equivalent to medium lesions (5–20 cm)
 - 5 mm device patency is equivalent to other sizes
 - Appears to have no dependence on device diameter in contrast to previous experience¹
- **Sizing is Critical**
 - Primary patency is significantly better when IFU sizing is not exceeded at the proximal edge
 - 88% versus 70% at 12 months ($p < .05$, sizing by Core Lab)
- **European VIASTAR Trial adds more Comparative Data**
 - Randomized trial of Bare Nitinol Stents versus the GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface for long SFA lesions

1. Saxon RR, Coffman JM, Gooding JM, Ponec DJ. Long-term patency and clinical outcome of the Viabahn Stent-Graft for femoropopliteal artery obstructions. *Journal of Vascular & Interventional Radiology* 2007;18(11):1341-1350.

Primary Patency in SFA Stenting

PRIMARY PATENCY IN SFA STENTING



Prospective multi-center studies with greater than two centers included. Trend line of randomized Bare Metal Stents (BMS's) only. Clinical study names, patient demographics, lesion characterization and patency definitions may differ among studies.

1. GORE® VIABAHN® Endoprosthesis
2. GORE® VIABAHN® with Heparin Bioactive Surface, labeled in some markets as the GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface
3. COOK® ZILVER® Stent
4. COOK® ZILVER® PTX® Stent

© 2015 W. L. Gore & Associates, Inc.

Bosiers M, Deloose K, Callaert J, et al. Results of the Protégé EverFlex 200-mm-long nitinol stent (ev3) in TASC C and D femoropopliteal lesions. *Journal of Vascular Surgery* 2011;54(4):1042-1050.



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface is known in some markets as the GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface.

Products listed may not be available in all markets.

Astron® Pulsar is a registered trademark of BIOTRONIK. SUPERA is part of The Leipzig Stent Registry. Zilver® PTX® is a registered trademark of Cook Medical. GORE®, PROPATEN, GORE-TEX®, PERFORMANCE BY DESIGN, VIABAHN®, and designs are trademarks of W. L. Gore & Associates. CARMEDA® and CBAS® are trademarks of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc. © 2011-2012, 2014-2015 W. L. Gore & Associates, Inc. AU0074-EU2 JULY 2015



VIPER
CLINICAL STUDY