

**GORE® EXCLUDER®**

Thoracoabdominal Branch  
Endoprosthesis

# BRIDGING THE GAP between thoracic and abdominal care

GET STARTED 



*Together, improving life*




**FIRST FDA APPROVED,**  
off-the-shelf endovascular  
solution for a range of  
complex visceral  
aortic cases.


# Meet the component you were looking for

Modular design and studied compatibility  
with other Gore commercial devices.


[WATCH DEPLOYMENT](#) 




**Aortic Component** positioned and  
partially deployed.



**Branch Components** delivered and  
deployed through pre-cannulated  
portals of Aortic Component.



**Aortic Component** deployment  
completed, followed by **Distal  
Bifurcated Component** and  
**Contralateral Leg Component**.



Seal zones ballooned  
as necessary.

[SEE THE DATA](#) 

# On-label solution for endovascular repair

- Thoracoabdominal aortic aneurysm patients with appropriate anatomy
- High-risk patients with pararenal aortic aneurysms and appropriate anatomy

Off-the-shelf availability for timely patient care

The benefits of minimally invasive care<sup>1,2</sup> without the time and effort to create custom grafts.



## Proven clinical performance, improved patient experience

TAMBE Pivotal Study through 12-month follow-up<sup>3</sup>

|       |                                                              |
|-------|--------------------------------------------------------------|
| 0%    | Lesion-related mortality                                     |
| 94%   | Freedom from aneurysm growth                                 |
| 92.6% | Freedom from failure of device effectiveness                 |
| 85.3% | Freedom from any device component occlusions                 |
| 100%  | Freedom from occlusion of aortic and iliac device components |
| 95.8% | Freedom from target visceral artery occlusion                |

### Real-world evidence (RWE)

Significant reduction in length of stay (LOS) for the endovascular procedure over open surgical repair.<sup>\*,4</sup>

15.8  
days

RWE average LOS for  
open surgical repair<sup>4</sup>

6.2  
days

RWE average LOS for  
endovascular repair<sup>4</sup>

4.9  
days ↓

Average LOS for  
TAMBE repair<sup>3</sup>



\* Data on file 2024; W. L. Gore & Associates, Inc; Flagstaff, AZ.

LEARN MORE



# Another first in branched technology

With dedicated support you can depend on in complex procedures.

- A deep reservoir of product knowledge.
- Our representatives support hundreds of aortic cases per year.
- Non-commissioned sales force – our focus is on outcomes.

---

Contact your Gore representative for more information.

---

## References

1. Bavaria JE, Appoo JJ, Makaroun MS, *et al*; Gore TAG Investigators. Endovascular stent grafting versus open surgical repair of descending thoracic aortic aneurysms in low-risk patients: a multicenter comparative trial. *Journal of Thoracic & Cardiovascular Surgery* 2007;133(2):369-377.
2. Arko FR 3rd, Murphy EH, Boyes C, *et al*. Current status of endovascular aneurysm repair: 20 years of learning. *Seminars in Vascular Surgery* 2012;25(3):131-135.
3. GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis [Instructions for Use]. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. MD193085.
4. PINC AI™ Applied Sciences, Premier Inc. PINC AI™ Healthcare Database: Data that informs and performs (White Paper). July 2023. <https://offers.pinc-ai.com/PINC-AI-Healthcare-Database-White-Paper-LP.html>



Consult Instructions  
for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

**INDICATIONS FOR USE IN THE U.S.:** The GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis is indicated for endovascular repair in patients with thoracoabdominal aortic aneurysms and high-surgical risk patients with pararenal aortic aneurysms who have appropriate anatomy as described below: Adequate iliac/femoral access and brachial/axillary access; Proximal (supraceliac) aortic neck treatment diameter range over 2 cm seal zone of 22–34 mm for aneurysms extending up to 6.5 cm or less above the origin of the most proximal branch vessel; Aortic neck angle  $\leq 60^\circ$  at the Aortic Component proximal seal zone; Iliac artery treatment diameter range of 8–25 mm and iliac artery seal zone length of at least 10 mm; Renal artery seal zone diameters between 4.0–10.0 mm; Celiac and superior mesenteric artery seal zone diameters between 5.0–12.0 mm;  $\geq 15$  mm seal zone length in renal arteries, superior mesenteric artery, and celiac artery; Visceral segment of aorta (3 cm proximal through 9.5 cm distal to the most proximal visceral artery) must be  $\geq 20$  mm in diameter. **CONTRAINDICATIONS:** The GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis is contraindicated in: Patients with known sensitivities or allergies to the TAMBE Device materials including ePTFE, FEP, nickel titanium alloy (Nitinol), stainless steel, and gold. Patients who have a condition that threatens to infect the graft. Patients with known hypersensitivity to heparin, including patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II and cannot receive the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis. <sup>Rx Only</sup>

Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available.

Products listed may not be available in all markets.

GORE, *Together, improving life*, EXCLUDER, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.  
© 2024 W. L. Gore & Associates, Inc. 241403251-EN MAY 2024

---

W. L. Gore & Associates, Inc.  
[goremedical.com](http://goremedical.com)

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673  
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

*Together, improving life*

