### **GORE® EXCLUDER®**


Thoracoabdominal Branch Endoprosthesis






## FIRST FDA-APPROVED, off-the-shelf endovascular solution for a range of complex visceral aortic cases.

**Aortic Component** positioned and partially deployed.



**Branch Components** delivered and deployed through pre-cannulated portals of Aortic Component.




**Aortic Component** deployment completed, followed by Distal **Bifurcated Component and** Contralateral Leg Component.

# Meet the component you were looking for

Modular design and studied compatibility with other Gore commercial devices.

WATCH DEPLOYMENT >





On-label clinically proven endovascular repair

 Thoracoabdominal aortic aneurysm patients with appropriate anatomy

 High-risk patients with pararenal aortic aneurysms and appropriate anatomy

Pivotal Study outcomes through 12 months<sup>1</sup>

0%	Lesion-related mortality
94%	Freedom from aneurysm growth
92.6%	Freedom from failure of device effectiveness
85.3%	Freedom from any device component occlusions
100%	Freedom from occlusion of aortic and iliac device components
95.8%	Freedom from target visceral artery occlusion

Off-the-shelf availability for timely patient care

The benefits of minimally invasive care<sup>2,3</sup> without the time and effort to create custom grafts.

4.9 days

Average length of stay for TAMBE repair<sup>1</sup>





# Another first in branched technology

With dedicated support you can depend on in complex procedures.

- A deep reservoir of product knowledge.
- Our representatives support hundreds of aortic cases per year.
- Non-commissioned sales force our focus is on outcomes.

Contact your Gore representative for more information.

#### References

1. GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis [Instructions for Use]. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. MD193085.

- 2. Bavaria JE, Appoo JJ, Makaroun MS, *et al.*; Gore TAG Investigators. Endovascular stent grafting versus open surgical repair of descending thoracic aortic aneurysms in low-risk patients: a multicenter comparative trial. *Journal of Thoracic & Cardiovascular Surgery* 2007;133(2):369-377.
- 3. Arko FR 3rd, Murphy EH, Boyes C, *et al*. Current status of endovascular aneurysm repair: 20 years of learning. *Seminars in Vascular Surgery* 2012;25(3):131-135.



INDICATIONS FOR USE IN THE U.S.: The GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis is indicated for endovascular repair in patients with thoracoabdominal aortic aneurysms and high-surgical risk patients with pararenal aortic aneurysms who have appropriate anatomy as described below: Adequate iliac/femoral access and brachial/axillary access; Proximal (supraceliac) aortic neck treatment diameter range over 2 cm seal zone of 22–34 mm for aneurysms extending up to 6.5 cm or less above the origin of the most proximal branch vessel; Aortic neck angle  $\leq$  60° at the Aortic Component proximal seal zone; Iliac artery treatment diameter range of 8–25 mm and iliac artery seal zone length of at least 10 mm; Renal artery seal zone diameters between 4.0–10.0 mm; Celiac and superior mesenteric artery seal zone diameters between 5.0–12.0 mm;  $\geq$  15 mm seal zone length in renal arteries, superior mesenteric artery, and celiac artery; Visceral segment of aorta (3 cm proximal through 9.5 cm distal to the most proximal visceral artery) must be  $\geq$  20 mm in diameter. **CONTRAINDICATIONS:** The GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis is contraindicated in: Patients with known sensitivities or allergies to the TAMBE Device materials including ePTFE, FEP, nickel titanium alloy (Nitinol), stainless steel, and gold. Patients who have a condition that threatens to infect the graft. Patients with known hypersensitivity to heparin, including patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II and cannot receive the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis.  $\Re_{\text{Cony}}$  Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available.

Products listed may not be available in all markets.

GORE, *Together, improving life*, EXCLUDER, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates. © 2024 W. L. Gore & Associates, Inc. 24AR1088-EN01 NOVEMBER 2024

W. L. Gore & Associates, Inc. goremedical.com

**Asia Pacific** +65 67332882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673 **United States** Flagstaff, AZ 86004 800 437 8181 928 779 2771

