

Small change. Big impact.

Heparin-coated shunts significantly reduce shunt-related mortality vs. non-heparin-coated shunts¹

Outcomes of a recent retrospective study comparing heparin-coated grafts with non-heparin-coated grafts utilized for conduits in systemic-to-pulmonary shunts in neonates and infants.¹

Outcomes	Non-heparin-coated (n = 73)	Heparin-coated (n = 69)	P-value
Mortality	14 (19.2%)	3 (4.3%)	0.006
Shunt-related mortality	12 (16.4%)	2 (2.9%)	0.007
Shunt occlusion / Thrombosis	12 (16.4%)	2 (2.9%)	0.007



GORE® PROPATEN® Vascular Graft configured for Pediatric Shunt

Features CBAS Heparin Surface, a coating demonstrated to provide lasting thromboresistance²⁻⁶



GORE® PROPATEN® Vascular Graft configured for Pediatric Shunt catalogue

Catalogue number	Internal diameter (mm)	Standard length (cm)
HPT030005A	3	5
HPT030010A	3	10
HPT030015A	3	15
HPT350005A	3.5	5
HPT350010A	3.5	10
HPT350015A	3.5	15
HPT040005A	4	5
HPT040010A	4	10
HPT040015A	4	15
HPT050005A	5	5
HPT050010A	5	10
HPT050015A	5	15
HPT060015A	6	15

INSURANCE

GORE® PROPATEN® Vascular Graft configured for Pediatric Shunt may be covered under similar coding as GORE-TEX® Vascular Graft configured for Pediatric Shunt. Consult your reimbursement resources for complete information.

Ask your Gore sales associate for inventory management of GORE-TEX® Vascular Graft configured for Pediatric Shunt



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

- Ashfaq A, Soroya MS, Iyengar A, Federman M, Reemtsen BL. Heparin-coated grafts reduce mortality in pediatric patients receiving systemic-to-pulmonary shunts. *Pediatric Cardiology* 2018;39(3):473-477.
- Gore S, Andersson J, Biran R, Underwood C, Riesenfeld J. Heparin surfaces: impact of immobilization chemistry on hemocompatibility and protein adsorption. *Journal of Biomedical Materials Research Part B: Applied Biomaterials* 2014;102(8):1817-1824.
- Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® BioActive Surface heparin immobilization. *European Journal of Vascular & Endovascular Surgery* 2003;25(5):432-437.
- Freeman J, Chen A, Weinberg RJ, Okada T, Chen C, Lin PH. Sustained thromboresistant bioactivity with reduced intimal hyperplasia of heparin-bonded PTFE Propaten Graft in a chronic canine femoral artery bypass model. *Annals of Vascular Surgery* 2018;49:295-303. <http://www.sciencedirect.com/science/article/pii/S0890509617310981>
- Biran R, Pond D. Heparin coatings for improving blood compatibility of medical devices. *Advanced Drug Delivery Reviews* 2017;112:12-23. <https://www.sciencedirect.com/science/article/pii/S0169409X16303210>
- Carmeda AB. Carmeda® BioActive Surface (also known as CBAS® Heparin Surface) Reference List. Upplands Väsby, Sweden: Carmeda AB; 2018. [Reference list]. MD166953

Indication for Use: GORE® PROPATEN® Vascular Grafts are intended for use as vascular prostheses for replacement or bypass of diseased vessels in patients suffering occlusive or aneurysmal diseases, in trauma patients requiring vascular replacement, for dialysis access, or for other vascular procedures. **CONTRAINDICATIONS** A. DO NOT use the GORE® PROPATEN® Vascular Graft in patients with known hypersensitivity to heparin, including those patients who have had a previous incidence of HIT type II. B. DO NOT use any configuration of GORE® PROPATEN® Vascular Grafts with Removable Rings, Non-Removable Rings or Integrated Rings for coronary artery bypass or cerebral reconstruction procedures. C. DO NOT use GORE® PROPATEN® Vascular Grafts as a patch. If cut and used as a patch, GORE® PROPATEN® Vascular Grafts may lack adequate transverse strength.

Important Information: This information is intended for education and awareness only. Patients should consult their physician for information on the risks associated with the devices and surgical procedures discussed in this website. All surgical procedures carry potential health risks. Not all patients will be candidates for treatment with these devices, and individual outcomes may vary. Always follow physician advice on your post-surgery care and recovery. Refer to the *Instructions for Use* for a complete list of warnings, precautions, contraindications, potential adverse events, and patient-specific information. Caution: U.S. law restricts use of this device on the order of a physician.

Products listed may not be available in all markets.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc. GORE®, GORE-TEX®, PROPATEN®, and designs are trademarks of W. L. Gore & Associates.

© 2019 W. L. Gore & Associates, Inc. AX1675-EN1 FEBRUARY 2019