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References

1. Moore J, Hegde S, El-Said J, et al; ACC IMPACT Steering Committee. Transcatheter device closure of atrial septal defects: a safety review. *JACC: Cardiovascular Interventions* 2013;6(5):433-442.
2. McElhinney DB, Quartermain MD, Kenny D, Alboliras E, Amin Z. Relative risk factors for cardiac erosion following transcatheter closure of atrial septal defects: a case-control study. *Circulation* 2016;133(18):1738-1746.
3. DiBardino DJ, McElhinney DB, Kaza AK, Mayer JE Jr. Analysis of the US Food and Drug administration manufacturer and user facility device experience database for adverse events involving Amplatzer Septal Occluder devices and comparison with the Society of Thoracic Surgery congenital cardiac surgery database. *Journal of Thoracic & Cardiovascular Surgery* 2009;137(6):1334-1341.
4. Thomson JDR, Qureshi SA. Device closure of secundum atrial septal defect's and the risk of cardiac erosion. *Echo Research & Practice* 2015;2(4): R73-R78.
5. Delaney JW, Li JS, Rhodes JF. Major complications associated with transcatheter atrial septal occluder implantation: a review of the medical literature and the manufacturer and user facility device experience (MAUDE) database. *Congenital Heart Disease* 2007;2(4):256-264.
6. El-Said HG, Moore JW. Erosion by the Amplatzer Septal Occluder: experienced operator opinions at odds with manufacturer recommendations?. *Catheterization & Cardiovascular Interventions* 2009;73(7):925-930.
7. Kitano M, Yazaki S, Sugiyama H, Ohtsuki S, Tomita H. Risk factors and predictors of cardiac erosion discovered from 12 Japanese patients who developed erosion after atrial septal defect closure using Amplatzer Septal Occluder. *Pediatric Cardiology* 2020;41(2):297-308.
8. Turner DR, Owada CY, Sang CJ Jr, Khan M, Lim DS. Closure of secundum atrial septal defects with the AMPLATZER Septal Occluder: a prospective, multicenter, post-approval study. *Circulation: Cardiovascular Interventions* 2017;10(8):e004212.
9. Taggart NW, Dearani JA, Hagler DJ. Late erosion of an Amplatzer Septal Occluder device 6 years after placement. *Journal of Thoracic & Cardiovascular Surgery* 2011;142(1):221-222.
10. Bergonti M, Toscano O, Teruzzi G, Trabattoni D. Never drop your guard down after atrial septal defect closure: a case report. *European Heart Journal. Case Reports* 2019;3(2):ytz094.
11. Amin Z, Hijazi ZM, Bass JL, Cheatham JP, Hellenbrand W, Kleinman CS. PFO closure complications from the AGA registry. *Catheterization & Cardiovascular Interventions* 2008;72(1):74-79.
12. Abe T, Tsukano S, Tosaka Y. Pericardial tamponade due to erosion of an Figulla Flex II device after closure of an atrial septal defect. *Catheterization & Cardiovascular Interventions* 2019;94(7):1003-1005.
13. Firouzi A, Siraj A, Masoud A, et al. Long-term outcome of atrial septal defect (ASD) repair by the Figulla Flex II ASD Occluder: a retrospective cohort study at a referral heart center in Iran. *Journal of Cardiology & Current Research* 2020;13(1):15-19.
14. Auriau J, Bouvaist H, Aaberge L, et al. Cardiac erosions after transcatheter atrial septal defect closure with the Figulla Flex device. *JACC: Cardiovascular Interventions* 2019;12(14):1397-1399.
15. Ikeda S, Imai M, Yoshizawa K, Miyamoto T, Fujiwara K, Sato Y. Cardiac erosion causing an aortic dissection after a catheter closure of an atrial septal defect: first late erosion case with the Figulla Flex II Septal Occluder. *Cardiovascular Intervention & Therapeutics* 2021;36(2):270-272.
16. Kister T, Dähnert I, Lurz P. Fatal erosion atrial septal defect device. *Catheterization & Cardiovascular Interventions* 2016;87(5):951-954.
17. Bier ML, Dhawan P, Shah SU, et al. Cardiac erosions with the Amplatzer Septal Occluder: adverse events in the Manufacturer and User Facility Device Experience (MAUDE) Database since the 2012 FDA review. *Structural Heart* 2021;5(1):85-89.
18. Crawford GB, Brindis RG, Krucoff MW, Mansalis BP, Carroll JD. Percutaneous atrial septal occluder devices and cardiac erosion: a review of the literature. *Catheterization & Cardiovascular Interventions* 2012;80(2):157-167.
19. Kumar P, Orford JL, Tobis JM. Two cases of pericardial tamponade due to nitinol wire fracture of a Gore Septal Occluder. *Catheterization & Cardiovascular Interventions* 2020;96(1):219-224.
20. W. L. Gore & Associates. *Safe and Effective Closure for Atrial Septal Defects with GORE® CARDIOFORM Septal Occluder and GORE® CARDIOFORM ASD Occluder*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2020. [Wire frame fracture analysis]. AZ1882-EN1.
21. W. L. Gore & Associates. *Safe and Effective Closure for Atrial Septal Defects. GORE® CARDIOFORM Septal Occluder, GORE® CARDIOFORM ASD Occluder*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2020. [Product brochure]. AZ0912-EN1.
22. Varotto L, Bregolin G, Pacanaro M, De Boni A, Bonanno C, Perini F. The closure of patent foramen ovale in preventing subsequent neurological events: a Bayesian network meta-analysis to identify the best device. *Cerebrovascular Diseases* 2020;49(2):124-134.
23. O'Byrne ML, Glatz AC, Sunderji S, et al. Prevalence of deficient retro-aortic rim and its effects on outcomes in device closure of atrial septal defects. *Pediatric Cardiology* 2014;35(7):1181-1190.

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© 2023 W. L. Gore & Associates, Inc. / 日本ゴア合同会社 22446850-JA FEBURARY 2023