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* Data on file 1999; W. L. Gore & Associates, Inc; Flagstaff, AZ.

References

1. Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® bioactive surface heparin immobilization. *European Journal of Vascular and Endovascular Surgery* 2003;25(5):432-437.
2. Biran R. Heparin Activity and Concentration Values from 52 and 96 Month Duration Clinical Explants of the GORE® PROPATEN® Vascular Graft. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2011. [Work plan]. WP104647.
3. Fisher JL, Thomson RC, Moore JW, Begovac PC. Functional parameters of thromboresistant heparinized e-PTFE vascular grafts. *Cardiovascular Pathology* 2002;11(1):42.
4. Heyligers JMM, Lisman T, Verhagen HJM, Weeterings C, de Groot PG, Moll FL. A heparin-bonded vascular graft generates no systemic effect on markers of hemostasis activation or detectable heparin-induced thrombocytopenia-associated antibodies in humans. *Journal of Vascular Surgery* 2008;47(2):324-329.
5. Lin PH, Chen C, Bush RL, Yao Q, Lumsden AB, Hanson SR. Small-caliber heparin-coated ePTFE grafts reduce platelet deposition and neointimal hyperplasia in a baboon model. *Journal of Vascular Surgery* 2004; 39(6):1322-1328.
6. Lin PH, Bush RL, Yao Q, Lumsden AB, Chen C. Evaluation of platelet deposition and neointimal hyperplasia of heparin-coated small-caliber ePTFE grafts in a canine femoral artery bypass model. *Journal of Surgical Research* 2004;118:45-52.
7. Kocsis JF, Llanos G, Holmer E. Heparin-coated stents. *Journal of Long-Term Effects of Medical Implants* 2000;10(1&2):19-45.
8. Chen C, Lumsden AB, Hanson SR. Local infusion of heparin reduces anastomotic neointimal hyperplasia in aortoiliac expanded polytetrafluoroethylene bypass grafts in baboons. *Journal of Vascular Surgery* 2000;31(2):354-363.
9. Heyligers JMM, Lisman T, Weeterings C, et al. Heparin immobilization reduces thrombogenicity on small-caliber ePTFE grafts. *Journal of Vascular Surgery* 2006;43:587-591.
10. Piffaretti G, Dorigo W, Castelli P, Pratesi C, Pulli R; PROPATEN Italian Registry Group. Results from a multicenter registry of heparin-bonded expanded polytetrafluoroethylene graft for above-the-knee femoropopliteal bypass. *Journal of Vascular Surgery* 2018;67(5):1463-1471.e1.
11. Samson RH, Morales R, Showalter DP, Lepore MR, Nair DG. Heparin-bonded expanded polytetrafluoroethylene femoropopliteal bypass grafts outperform expanded polytetrafluoroethylene grafts without heparin in a long-term comparison. *Journal of Vascular Surgery* 2016;64(3):638-647.
12. Lindholt JS, Gottschalksen B, Johannesen N, et al. The Scandinavian Propaten® Trial – 1-year patency of PTFE vascular prostheses with heparin-bonded luminal surfaces compared to ordinary pure PTFE vascular prostheses – a randomised clinical controlled multi-centre trial. *European Journal of Vascular & Endovascular Surgery* 2011;41(5):668-673.
13. Daenens K, Schepers S, Fourneau I, Houthoofd S, Nevelsteen A. Heparin-bonded ePTFE grafts compared with vein grafts in femoropopliteal and femorocrural bypasses: 1- and 2-year results. *Journal of Vascular Surgery* 2009;49(5):1210-1216.
14. Peeters P, Verbist J, Deloose K, Bosiers M. Will heparin-bonded PTFE replace autologous venous conduits in infrapopliteal bypass? *Italian Journal of Vascular & Endovascular Surgery* 2008;15(3):143-148.
15. Bosiers M, Deloose K, Verbist J, et al. Heparin-bonded expanded polytetrafluoroethylene vascular graft for femoropopliteal and femorocrural bypass grafting: 1-year results. *Journal of Vascular Surgery* 2006;43(2):313-319.
16. Kaisar J, Chen A, Cheung M, Kfoury E, Bechara C, Lin P. Comparison of propaten heparin-bonded vascular graft with distal anastomotic patch versus autogenous saphenous vein graft in tibial artery bypass. *Vascular* 2018;26(2):117-125.
17. Uhl C, Hock C, Betz T, Greindl M, Töpel I, Steinbauer M. Comparison of venous and HePTFE tibial and peroneal bypasses in critical limb ischemia patients unsuitable for endovascular revascularization. *Vascular* 2015;23(6):607-613.
18. Neville RF, Babrowicz J, Amdur R, Sidawy A. Heparin-bonded ePTFE (Propaten): is it as good as autologous vein for tibial bypass? *Italian Journal of Vascular & Endovascular Surgery* 2014;21(1)49-56.
19. Monaca V, Battaglia G, Turiano SA, Tringale R, Catalfamo S. Subpopliteal revascularization. Criteria analysis for use of E-PTFE (Propaten®) as first choice conduit. *Italian Journal of Vascular & Endovascular Surgery* 2013;20(3):165-169.
20. Dorigo W, Raffaele P, Piffaretti G. Results from an Italian multicentric registry comparing heparin-bonded ePTFE graft and autologous saphenous vein in below-knee femoro-popliteal bypasses. *Journal of Cardiovascular Surgery* 2012;53(2):187-193.
21. Pulli R, Dorigo W, Castelli P, et al. Propaten Italian Registry Group. Midterm results from a multicenter registry on the treatment of infrainguinal critical limb ischemia using a heparin-bonded ePTFE graft. *Journal of Vascular Surgery* 2010;51(5):1167-1177.e1.
22. Hugl B, Nevelsteen A, Daenens K, et al. PEPE II Study Group. PEPE II – a multicenter study with an end-point heparin-bonded expanded polytetrafluoroethylene vascular graft for above and below knee bypass surgery: determinants of patency. *Journal of Cardiovascular Surgery* 2009;50(2):195-203.

References

23. Lösel-Sadée H, Alefelder C. Heparin-bonded expanded polytetrafluoroethylene graft for infragenicular bypass: five-year results. *Journal of Cardiovascular Surgery* 2009;50(3):339-343.
24. Belch JJ, Dormandy J, Biasi GM, *et al*; CASPAR Writing Committee. Results of the randomized, placebo-controlled clopidogrel and acetylsalicylic acid in bypass surgery for peripheral arterial disease (CASPAR) trial. *Journal of Vascular Surgery* 2010;52(4):825-833.
25. Warkentin TE. An overview of the heparin-induced thrombocytopenia syndrome. *Seminars in Thrombosis & Hemostasis* 2004;30(3):273-283.
26. Franchini M. Heparin-induced thrombocytopenia: an update. *Thrombosis Journal* 2005;3(1):14-18.
27. Jang IK, Hursting MJ. When heparins promote thrombosis: review of heparin-induced thrombocytopenia. *Circulation* 2005;111(20):2671-2683.
28. Chong BH. Heparin-induced thrombocytopenia. *Journal of Thrombosis & Haemostasis* 2003;1(7):1471-1478.
29. Lubenow N, Kempf R, Eichner A, Eichler P, Carlsson LE, Greinacher A. Heparin-induced thrombocytopenia: temporal pattern of thrombocytopenia in relation to initial use or reexposure to heparin. *Chest* 2002;122(1):37-42.
30. Warkentin TE, Kelton JG. Temporal aspects of heparin-induced thrombocytopenia. *New England Journal of Medicine* 2001;344(17):1286-1292.
31. Warkentin TE, Griener A, Koster A, Lincoff AM. Treatment and prevention of heparin-induced thrombocytopenia: American College of Chest Physicians evidence-based clinical practice guidelines (8th edition). *Chest* 2008;133:3405-3805.
32. Jackson MR, Neilson WJ, Lary M, Baay P, Web K, Clagett GP. Delayed-onset heparin-induced thrombocytopenia and thrombosis after intraoperative heparin anticoagulation-four case reports. *Vascular & Endovascular Surgery* 2006;40(1):67-70.
33. Rice L, Attisha WK, Drexler A, Francis JL. Delayed-onset heparin-induced thrombocytopenia. *Annals of Internal Medicine* 2002;136(3):210-215.
34. Smythe MA, Stephens JL, Mattson JC. Delayed-onset heparin-induced thrombocytopenia. *Annals of Emergency Medicine* 2005;45(4):417-419.
35. Warkentin TE, Kelton JG. Delayed-onset heparin-induced thrombocytopenia and thrombosis. *Annals of Internal Medicine* 2001;135(7):502-506.
36. Kasirajan K. Outcomes following heparin-induced thrombocytopenia in patients with heparin bonded vascular grafts. Abstract presented at the 2011 Vascular Annual Meeting; June 16-18, 2011; Chicago, IL. *Journal of Vascular Surgery* 2011;53(6)Supplement 1:17S-18S. Abstract SS4.

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© 2025 W. L. Gore & Associates, Inc. / 日本ゴア合同会社 無断転載・複製を禁じます。記載の商標は、Goreグループの関連会社または各権利者の商標です。