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| 4th Step: CDT               |                        |
|-----------------------------|------------------------|
| Continuous                  | Intermittent           |
| - Saline 1000ml             | - Urokinase 6000 units |
| - Heparin 10 - 16 ml        | 4 times/day            |
| - Argatroban 20mg           | - Cefazolin 1g         |
| 40ml/hr                     | 2 times/day            |
| <b>Confirm L/D everyday</b> |                        |
| <b>Periods: 1-2 days</b>    |                        |

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## Reference

- Banerjee S, Sarode K, Mohammad A, et al. Femoropopliteal artery stent thrombosis: report from the Excellence in Peripheral Artery Disease Registry. *Circulation: Cardiovascular Interventions* 2016;9(2):e002730.
- Iida O, Takahara M, Soga Y, et al; VANQUISH Investigators. One-year outcomes of heparin-bonded stent-graft therapy for real-world femoropopliteal lesions and the association of patency with the prothrombotic state based on the prospective, observational, multicenter Viabahn Stent-Graft Placement for Femoropopliteal Diseases Requiring Endovascular Therapy (VANQUISH) Study. *Journal of Endovascular Therapy* 2021;28(1):123-131.
- Saxon RR, Chervu A, Jones PA, et al. Heparin-bonded, expanded polytetrafluoroethylene-lined stent graft in the treatment of femoropopliteal artery disease: 1-year results of the VIPER (Viabahn Endoprosthesis with Heparin Bioactive Surface in the Treatment of Superficial Femoral Artery Obstructive Disease) trial. *Journal of Vascular & Interventional Radiology* 2013;24(2):165-173.
- Ohki T, Kichikawa K, Yokoi H, et al. Outcomes of the Japanese multicenter Viabahn trial of endovascular stent grafting for superficial femoral artery lesions. *Journal of Vascular Surgery* 2017;66(1):130-142.e1.
- Kurata N, Iida O, Asai M, et al. Factors influencing in-stent occlusion after femoropopliteal artery stent placement with intravascular ultrasound evaluation. *Journal of Vascular & Interventional Radiology* 2020;31(2):213-220-224.
- Patel SR, Hughes CO, Jones KG, et al. A systematic review and meta-analysis of endovascular popliteal aneurysm repair using the Hemobahn/Viabahn Stent-Graft. *Journal of Endovascular Therapy* 2015;22(3):330-337.
- Tomoi Y, Soga Y, Iida O, Fujihara M, Ando K. Impact of Cilostazol administration on femoropopliteal in-stent restenosis. *Journal of Endovascular Therapy* 2017;24(5):640-646.
- Soga Y, Hamasaki T, Eda Hiro R, et al; STOP-IC Investigators. Sustained effectiveness of Cilostazol after endovascular treatment of femoropopliteal lesions: midterm follow-up from the sufficient treatment of peripheral intervention by Cilostazol (STOP-IC) study. *Journal of Endovascular Therapy* 2018;25(3):306-312.
- Geng DF, Liu M, Jin DM, Wu W, Deng J, Wang JF. Cilostazol-based triple antiplatelet therapy compared to dual antiplatelet therapy in patients with coronary stent implantation: a meta-analysis of 5,821 patients. *Cardiology* 2012;122(3):148-157.

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### Reference

1. Iida O, Takahara M, Soga Y, et al; VANQUISH Investigators. One-year outcomes of heparin-bonded stent-graft therapy for real-world femoropopliteal lesions and the association of patency with the prothrombotic state based on the prospective, observational, multicenter Viabahn Stent-Graft Placement for Femoropopliteal Diseases Requiring Endovascular Therapy (VANQUISH) Study. *Journal of Endovascular Therapy* 2021;28(1):123-131.

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© 2022 W. L. Gore & Associates, Inc. / 日本ゴア合同会社 21217373-JA APRIL 2022
W. L. Gore & Associates, Inc.