

**GORE® TAG®**

Conformable Thoracic Stent Graft  
with ACTIVE CONTROL System

# WHERE CONFORMABILITY MEETS CONTROL

Combining the precise control of an intuitive delivery system with the exceptional conformability of the GORE® TAG® Conformable Thoracic Stent Graft, clinicians can achieve predictable patient outcomes.<sup>1-3</sup>

- Staged deployment enables refinement of placement
- Optional angulation control at intermediate and full diameter
- Hemodynamic stability throughout the deployment
- Expanded large-diameter taper sizes



*Together, improving life*




## GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System

Catalogue number	Intended aortic diameter (mm)	Proximal diameter (mm)	Distal diameter (mm)	Endoprosthesis length (cm)	Recommended GORE® DRYSEAL Flex Introducer Sheath size (Fr)	Oversizing range (%)	Partially uncovered stent length (mm)
<b>Straight Devices</b>							
TGM212110	16–19.5	21	21	10	18	8–31	3
TGM262610	19.5–24	26	26	10	20	8–33	4
TGM282810	22–26	28	28	10	20	8–27	4
TGM282815	22–26	28	28	15	20	8–27	4
TGMR313110	24–29	31	31	10	20	7–29	4
TGMR313115	24–29	31	31	15	20	7–29	4
TGMR313120	24–29	31	31	20	20	7–29	4
TGM343410	27–32	34	34	10	22	6–26	5
TGM343415	27–32	34	34	15	22	6–26	5
TGM343420	27–32	34	34	20	22	6–26	5
TGMR373710	29–34	37	37	10	22	9–28	5
TGMR373715	29–34	37	37	15	22	9–28	5
TGMR373720	29–34	37	37	20	22	9–28	5
TGMR404010	31–37	40	40	10	22	8–29	6
TGMR404015	31–37	40	40	15	22	8–29	6
TGMR404020	31–37	40	40	20	22	8–29	6
TGM454510	34–42	45	45	10	24	7–32	6.5
TGM454515	34–42	45	45	15	24	7–32	6.5
TGM454520	34–42	45	45	20	24	7–32	6.5
<b>Taper Devices*</b>							
TGM262110	19.5–24/16–19.5	26	21	10	20	8–33	4
TGMR312610	24–29/19.5–24	31	26	10	20	7–33	4
TGM342815	27–32/22–26	34	28	15	22	6–27	5
TGMR373115	29–34/24–29	37	31	15	22	7–29	5
TGMR403415	31–37/27–32	40	34	15	22	6–29	6
TGM453715	34–42/29–34	45	37	15	24	7–32	6.5

\* Some sizes not available in all markets.

1. OP869 Farber MA, Krishnasastry KV, Desai N, Starnes BW, Matsumura JS, Tohill BC; TAG 08-02 Clinical Trial Investigators. Five-year outcomes with Conformable GORE® TAG® endoprosthesis used in traumatic aortic transections. *Annals of Thoracic Surgery*. In press. <https://www.sciencedirect.com/science/article/pii/S0003497521010134>
2. OP871 Jordan WD, Desai N, Letter AJ, Matsumura JS. Long-term outcomes of the CTAG thoracic endograft in a prospective multicenter trial. *Journal of Vascular Surgery* 2021;74(5):1491-1498. <https://www.sciencedirect.com/science/article/pii/S0741521421007485>
3. OP1025 Desai N, Tohill B, Matsumura J, Cambria RP. A multicenter clinical trial on the five-year outcomes following treatment of acute, complicated type B aortic dissection with a conformable stent graft. *Annals of Vascular Surgery* 2025;113:1-12. A Multicenter Clinical Trial on the Five-Year Outcomes Following Treatment of Acute, Complicated Type B Aortic Dissection With a Conformable Stent Graft - ScienceDirect.

 Consult Instructions for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

**INDICATIONS FOR USE IN THE U.S.:** The GORE® TAG® Conformable Thoracic Stent Graft is intended for endovascular repair of all lesions of the descending thoracic aorta, including: Isolated lesions in patients who have appropriate anatomy, including: adequate iliac/femoral access, aortic inner diameter in the range of 16–42 mm, ≥ 20 mm non-aneurysmal aorta proximal and distal to the lesion; Type B dissections in patients who have appropriate anatomy, including: adequate iliac/femoral access, ≥ 20 mm landing zone proximal to the primary entry tear; proximal extent of the landing zone must not be dissected, diameter at proximal extent of proximal landing zone in the range of 16–42 mm. **CONTRAINDICATIONS:** Patients with known sensitivities or allergies to the device materials; patients who have a condition that threatens to infect the graft. Refer to *Instructions for Use* at [eifu.goremedical.com](http://eifu.goremedical.com) for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available.  only  
 Products listed may not be available in all markets.

GORE, *Together, improving life*, ACTIVE CONTROL, DRYSEAL, TAG and designs are trademarks of W. L. Gore & Associates.  
 © 2025 W. L. Gore & Associates, Inc. 25AR1004-EN01 MARCH 2025

W. L. Gore & Associates, Inc.  
[goremedical.com](http://goremedical.com)

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673  
 United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

