


Reline

with Confidence.




PERFORMANCE through innovation

THE CONTINUING EVOLUTION OF A REVOLUTIONARY DEVICE



VIABAHN®
ENDOPROSTHESIS

HEPARIN
BIOACTIVE SURFACE



Gore continues
to evolve the

GORE® VIABAHN®
Endoprosthesis,

demonstrating our

commitment to

providing our

customers with


innovative products.




▶ A New Beginning for the Proximal End

- New precision laser trimming technology enables manufacturing change
- Excess material at the proximal edge removed
- Contoured edge may improve flow dynamics at proximal end
- The *Instructions for Use*, including sizing and placement recommendations, remain unchanged

Canine In Vivo IVUS Examples



Non-contoured edge



Contoured edge

Excess material removed at the device margin of the contoured edge compared to a non-contoured edge

▶ The Endoluminal SFA Bypass

Cover with Confidence

Covers and excludes the diseased irregular tissue of the arterial wall

ePTFE Lining

Provides barrier to in-stent restenosis

Nitinol Stent

Conformable yet durable

Heparin-Bonded Surface

Intended to provide sustained thromboresistance

CBAS® Heparin Surface

- Intended to provide a thromboresistant surface
- Sustained bioactivity*
- Proprietary end-point covalent bonding

Gore® VIABAHN® Endoprosthesis with Heparin Bioactive Surface



The bioactive luminal surface of a 5 mm diameter Gore® VIABAHN® Endoprosthesis with Heparin Bioactive Surface appears free of thrombus after two hours in an in vitro blood loop model.

Control Endoprosthesis




The non-bioactive luminal surface of a control endoprosthesis (5 mm diameter) appears covered with thrombus after 90 minutes in the same blood loop model (data on file).




Smaller Profiles

Reduced delivery profile for 5 – 8 mm devices


DIAMETER (mm)	PROFILE (Fr)	GUIDEWIRE
5	6	0.014" or 0.018"
6	6	0.014" or 0.018"
7	7	0.014" or 0.018"
8	7 ⁴	0.014" or 0.018"



New catheter design maintains pushability and trackability of previous generation




► Sustained Bioactivity




Postoperative Week

Long-term Heparin Activity of Explanted Heparin-bonded ePTFE Vascular Grafts in a Canine Model*

* Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® bioactive surface heparin immobilization. *European Journal of Vascular and Endovascular Surgery* 2003;25(5):432-437.



Proprietary Covalent End-point Bonding



Covalent end-point bonding allows the heparin to extend into the bloodstream, keeping the active site bioavailable, unlike a non-permanent bond that can be washed away in the bloodstream.

- The anticoagulant function of heparin is dependent on the bioavailability of an active site within the molecule.
- Some methods of covalent heparin bonding damage and / or obstruct the active site, and hence destroy the heparin's anticoagulant activity.
- The CBAS® Heparin Surface of the GORE® PROPATEN® Vascular Graft consists of a proprietary covalent end-point bond that preserves the active site, thus retaining heparin's anticoagulant activity.

Sizing Table

0.035" Guidewire Compatibility

Device Sizing		Introducer Sheath (Fr)				
ENDOPROSTHESIS LABELED DIAMETER ¹ (mm)	RECOMMENDED VESSEL DIAMETER ² (mm)	2.5 cm DEVICE LENGTH ¹	5 cm DEVICE LENGTH ¹	10 cm DEVICE LENGTH ¹	15 cm DEVICE LENGTH ¹	RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP ³ (mm)
5	4.0 – 4.7	7	7	7	7	5
6	4.8 – 5.5	7	7	7	7	6
7	5.6 – 6.5	8	8	8	8	7
8	6.6 – 7.5	8	8	8	8	8
9	7.6 – 8.5	–	9	9	9	9
10	8.6 – 9.5	11 ⁴	11 ⁴	11 ⁴	11 ⁴	10
11	9.6 – 10.5	11	11	11	–	12
13	10.6 – 12.0	12	12	12	–	14

0.018" Guidewire Compatibility

Device Sizing		Introducer Sheath (Fr)				
ENDOPROSTHESIS LABELED DIAMETER ¹ (mm)	RECOMMENDED VESSEL DIAMETER ² (mm)	2.5 cm DEVICE LENGTH ¹	5 cm DEVICE LENGTH ¹	10 cm DEVICE LENGTH ¹	15 cm DEVICE LENGTH ¹	RECOMMENDED BALLOON DIAMETER FOR DEVICE TOUCH-UP ³ (mm)
5	4.0 – 4.7	6	6	6	6	5
6	4.8 – 5.5	6	6	6	6	6
7	5.6 – 6.5	7	7	7	7	7
8	6.6 – 7.5	7	7	7	7	8

¹ Labeled device diameters and lengths are nominal.

² Recommended endoprosthesis compression within the vessel is approximately 5 – 20%.

³ For the 11 and 13 mm diameter devices, balloon inflation pressure should not exceed 8 atm.

⁴ The 10 mm diameter device is compatible with the following 10 Fr introducer sheaths: Cordis AVANTI® Sheath Introducer, Boston Scientific SUPER SHEATH Introducer Sheath, B. Braun INTRADYN Tear-Away Introducer Sheath. The 8 mm x 25 cm device is not compatible with the 7 Fr COOK® CHECK-FLO® FLEXOR® Sheath.



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com

* Begovac PC, Thomson RC, Fisher JL, Hughson A, Gällhagen A. Improvements in GORE-TEX® Vascular Graft performance by Carmeda® bioactive surface heparin immobilization. *European Journal of Vascular and Endovascular Surgery* 2003;25(5):432-437.

Products listed may not be available in all markets.

SUPER SHEATH is a trademark of Boston Scientific Corporation. INTRADYN is a trademark of B. Braun Interventional Systems Inc. AVANTI® is a trademark of Cordis Corporation. COOK®, CHECK-FLO®, and FLEXOR® are trademarks of Cook Medical, Inc. GORE®, GORE-TEX®, HEMOBAHN®, PERFORMANCE THROUGH INNOVATION, TIP to HUB, VIABAHN®, and designs are trademarks of W. L. Gore & Associates. CARMEDA® and CBAS® are trademarks of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates. © 2015 W. L. Gore & Associates, Inc. AU0702-EN1 NOVEMBER 2015