



Revolutionary device

Proven performance

STUDY DESIGN

- Prospective, multi-center
- 120 limbs at 12 sites
- Lesion length ≥ 5 cm
- Primary patency by duplex (PSVR < 2.5)
- Independent Core Lab vessel sizing

LESION CHARACTERISTICS

- Average lesion length 19 cm
- 56% Chronic Total Occlusion

RESULTS

- One-year primary patency 73 percent (Figure 1)
- Patency independent of device diameter (5 mm device same as 6 mm device) and lesion length
- Primary patency is improved ($p < 0.05$) when IFU sizing not exceeded (Figure 2)

* Saxon RR, Chervu A, Jones PA, et al. Heparin-bonded, expanded polytetrafluoroethylene-lined stent graft in the treatment of femoropopliteal artery disease: 1-year results of the VIPER (Viabahn Endoprosthesis with Heparin Bioactive Surface in the Treatment of Superficial Femoral Artery Obstructive Disease) Trial. *Journal of Vascular & Interventional Radiology* 2013;24(2):165-173.

W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.

GORE®, VIABAHN®, and designs are trademarks of W. L. Gore & Associates.
© 2017 W. L. Gore & Associates, Inc. AW1470-EN1 NOVEMBER 2017

ONE-YEAR RESULTS FOR THE GORE VIPER CLINICAL STUDY*

Principle Investigator: Richard Saxon, MD

The first controlled, multi-center, prospective study evaluating GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface and Contoured Proximal Edge




Figure 1. One-year device patency




Figure 2. Primary patency versus oversizing at the proximal edge



VIPER

CLINICAL STUDY