Protesi preperitoneale
GORE® SYNECOR

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Bibliografia

1. W. L. Gore & Associates, Inc. *Clinical Evaluation Report for GORE® SYNECOR Preperitoneal Biomaterial*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Relazione sulla valutazione clinica – CER]. MD188925. Rev 1.
2. Gokcal F, Omar Yusef Kudsi OY. Are morbid obese (class–III) patients at high risk for postoperative complications after robotic ventral hernia repair? A propensity score matching analysis. *Laparoscopic Endoscopic Surgical Science* 2020;27(2):37-47.
3. Kudsi OY, Chang K, Bou-Ayash N, Gokcal F. Hybrid robotic hernia repair for incisional hernias: perioperative and patient-reported outcomes. *Journal of Laparoendoscopic & Advanced Surgical Techniques* 2021;31(5):570-578.
4. Kudsi OY, Kaoukabani G, Bou-Ayash N, et al. *Quality of life and surgical outcomes of robotic retromuscular ventral hernia repair using a new hybrid mesh reinforcement*. *Hernia* 2022;26(3):881-888.
5. Rios-Diaz A, Hitchner M, Christopher AN, Broach R, Cunning JR, Fischer JP. Early clinical and patient-reported outcomes of a new hybrid mesh for incisional hernia repair. *Journal of Surgical Research* 2021;265:49-59.
6. Petro CC, Nahabet EH, Criss CN, et al. Central failures of lightweight monofilament polyester mesh causing hernia recurrence: a cautionary note. *Hernia* 2015;19(1):155-159.
7. Cobb WS, Warren JA, Ewing JA, Burnikel A, Merchant M, Carbonell AM. Open retromuscular mesh repair of complex incisional hernia: predictors of wound events and recurrence. *Journal of the American College of Surgeons* 2015;220(4):606-613.
8. Warren JA, McGrath SP, Hale AL, Ewing JA, Carbonell AM 2nd, Cobb WS 4th. Patterns of recurrence and mechanisms of failure after open ventral hernia repair with mesh. *American Surgeon* 2017;83(11):1275-1282.
9. W. L. Gore & Associates, Inc. *GORE SYNECOR Preperitoneal Biomaterial Design Control (DC) Matrix*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Matrice di controllo della progettazione – Matrice DC]. MD187124. Rev 2.
10. Ventral Hernia Working Group; Breuing K, Butler CE, et al. Incisional ventral hernias: review of the literature and recommendations regarding the grading and technique of repair. *Surgery* 2010;148(3):544-558.
11. Crawford N. *Assessment of Vascularity via Micro CT in Various Patch Devices*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Rapporto finale dello studio]. 2344TL.
12. Berman A. *Evaluation of Plexus with no film and ETHICON PHYSIOMESH® in a 30-day rabbit subcutaneous model*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2015. [Protocollo dello studio]. 2336SC.
13. W. L. Gore & Associates, Inc. *Plexus Knit PQ Validation Report*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2022. [Rapporto di convalida]. MD145325. Rev 5.
14. Olson TB. *Competitive Hernia Device Strength Evaluation*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Piano di lavoro]. WP108484.
15. Sharkawy AA, Klitzman B, Truskey GA, Reichert WM. Engineering the tissue which encapsulates subcutaneous implants. II. Plasma-tissue exchange properties. *Journal of Biomedical Materials Research* 1998;40(4):586-597.
16. Holt DJ, Grainger DW. Host response to biomaterials. In: Hollinger JO, ed. *An Introduction to Biomaterials*. 2nd ed. Boca Raton, FL: CRC Press; 2012;6:91-118.
17. Berman A. *Evaluation of Plexus without film and ETHICON PHYSIOMESH® in a 180-day rabbit subcutaneous model*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2015 [Protocollo dello studio]. 2338SC.
18. Berman A. *A Rabbit Model for the Biomechanical and Histological Assessment of Suture Line Wound Healing*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2008. [Protocollo dello studio]. 1978SC.
19. W. L. Gore & Associates, Inc. *GORE® BIO-A® Tissue Reinforcement Design Control (DC) Matrix*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Matrice di controllo della progettazione – Matrice DC]. MD180587. Rev 7.
20. Morales-Conde S, Flores M, Fernandez V, Morales-Mendez S. *Bioabsorbable vs polypropylene plug for the “Mesh and Plug” inguinal hernia repair*. Poster presentato al 9th Annual Meeting of the American Hernia Society. 9-12 febbraio 2005; San Diego, CA.
21. LeBlanc KA, Heniford BT, Voeller GR. Innovations in ventral hernia repair. Materials and techniques to reduce MRSA and other infections. *Contemporary Surgery* 2006;62(4)Supplement:1-8.
22. Clinger LR. *PTFE Knit Microbial Placement*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2018. [Piano di lavoro]. WP110158.
23. Lake SP, Ray S, Zihni AM, Thompson DM Jr, Gluckstein J, Deeken CR. Pore size and pore shape – but not mesh density – alter the mechanical strength of tissue ingrowth and host tissue response to synthetic mesh materials in a porcine model of ventral hernia repair. *Journal of the Mechanical Behavior of Biomedical Materials* 2015;42:186-197.
24. Klinge U, Klosterhalfen B, Birkenhauer V, Junge K, Conze J, Schumpelick V. Impact of polymer pore size on the interface scar formation in a rat model. *Journal of Surgical Research* 2002;103(2):208-214.
25. Zhu LM, Schuster P, Klinge U. Mesh implants: an overview of crucial mesh parameters. *World Journal of Gastrointestinal Surgery* 2015;7(10):226-236.
26. Klinge U, Klosterhalfen B, Conze J, et al. Modified mesh for hernia repair that is adapted to the physiology of the abdominal wall. *European Journal of Surgery* 1998;164(12):951-960.
27. Socea B, Socea LI, Bratu OG, et al. Recurrence rates and mesh shrinkage after polypropylene vs. polyester mesh hernia repair in complicated hernias. *Materiale Plastice* 2018;55(1):79-81.
28. Klosterhalfen B, Junge K, Klinge U. The lightweight and large porous mesh concept for hernia repair. *Expert Reviews in Medical Devices* 2005;2(1):103-117.
29. Jenkins ED, Yip M, Melman L, Frisella MM, Matthews BD. Informed consent: cultural and religious issues associated with the use of allogeneic and xenogeneic mesh products. *Journal of the American College of Surgeons* 2010;210(4):402-410.
30. Lighter M, Roberts J. The role of biosynthetic mesh in abdominal wall hernia repair in the setting of obesity, recurrence and high risk patients. Presentato al 2019 AHS Annual Meeting; 11-14 marzo 2019; Las Vegas, NV. *Hernia* 2019;23 (Supplement 1): S96-S97. P-1256.
31. Landry M, Grimesly L, Gorman B, et al. Clinical Quality Improvement (CQI): evaluating the real work use of a new hybrid hernia mesh. Presentato al 2018 International Hernia Congress Program; 12 marzo, 2018; Miami Beach, FL.
32. Kudsi OY, Gokcal F. Lateral approach totally extraperitoneal (TEP) robotic retromuscular ventral hernia repair. *Hernia* 2021;25(1):211-222.
33. W. L. Gore & Associates. *Assessment of GORE® SYNECOR Biomaterial in Focused Patient Populations and Long-Term Application*. Identificatore NLM: NCT05094089. Data di pubblicazione: 26 ottobre 2021. Data di aggiornamento: 20 aprile 2022. Data di accesso: 19 maggio 2022. <https://clinicaltrials.gov/ct2/show/NCT05094089>
34. Yang XF, Liu JL. Anatomy essentials for laparoscopic inguinal hernia repair. *Annals of Translational Medicine* 2016;4(19):372.
35. Claus, C., Furtado, M., Malcher, F. et al. Ten golden rules for a safe MIS inguinal hernia repair using a new anatomical concept as a guide. *Surgical Endoscopy* 2020;34:1458-1464.

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© 2024 W. L. Gore & Associates GmbH 4GP2037-IT01 SETTEMBRE 2024
OPTILENE è un marchio commerciale di Aesculap AG.
BD, BARD e 3DMAX sono marchi commerciali di C.R. Bard, Inc.
TELA BIO e OVITEX sono marchi commerciali di TELA Bio, Inc.