GORE® SYNECOR

Präperitoneales Biomaterial

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Literatur

1. W. L. Gore & Associates, Inc. *Clinical Evaluation Report for GORE® SYNECOR Preperitoneal Biomaterial*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Klinischer Bewertungsbericht (Clinical Evaluation Report, CER)]. MD188925. 1. Fassung.
2. Gokcal F, Omar Yusef Kudsi OY. Are morbid obese (class–III) patients at high risk for postoperative complications after robotic ventral hernia repair? A propensity score matching analysis. *Laparoscopic Endoscopic Surgical Science* 2020;27(2):37–47.
3. Kudsi OY, Chang K, Bou-Ayash N, Gokcal F. Hybrid robotic hernia repair for incisional hernias: perioperative and patient-reported outcomes. *Journal of Laparoendoscopic & Advanced Surgical Techniques* 2021;31(5):570–578.
4. Kudsi OY, Kaoukabani G, Bou-Ayash N, et al. *Quality of life and surgical outcomes of robotic retromuscular ventral hernia repair using a new hybrid mesh reinforcement*. *Hernia* 2022;26(3):881–888.
5. Rios-Diaz A, Hitchner M, Christopher AN, Broach R, Cuning JR, Fischer JP. Early clinical and patient-reported outcomes of a new hybrid mesh for incisional hernia repair. *Journal of Surgical Research* 2021;265:49–59.
6. Petro CC, Nahabet EH, Criss CN, et al. Central failures of lightweight monofilament polyester mesh causing hernia recurrence: a cautionary note. *Hernia* 2015;19(1):155–159.
7. Cobb WS, Warren JA, Ewing JA, Burnikel A, Merchant M, Carbonell AM. Open retromuscular mesh repair of complex incisional hernia: predictors of wound events and recurrence. *Journal of the American College of Surgeons* 2015;220(4):606–613.
8. Warren JA, McGrath SP, Hale AL, Ewing JA, Carbonell AM 2nd, Cobb WS 4th. Patterns of recurrence and mechanisms of failure after open ventral hernia repair with mesh. *Surgeon* 2017;83(11):1275–1282.
9. W. L. Gore & Associates, Inc. *GORE SYNECOR Preperitoneal Biomaterial Design Control (DC) Matrix*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Design Control Matrix – DC Matrix]. MD187124. 2. Fassung.
10. Ventral Hernia Working Group; Breuing K, Butler CE, et al. Incisional ventral hernias: review of the literature and recommendations regarding the grading and technique of repair. *Surgery* 2010;148(3):544–558.
11. Crawford N. *Assessment of Vascularity via Micro CT in Various Patch Devices*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Abschließender Studienbericht]. 2344TL.
12. Berman A. *Evaluation of Plexus with no film and ETHICON PHYSIOMESH® in a 30-day rabbit subcutaneous model*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2015. [Studienprotokoll]. 23365C.
13. W. L. Gore & Associates, Inc. *Plexus Knit PQ Validation Report*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2022. [Validierungsbericht]. MD145325. 5. Fassung.
14. Olson TB. *Competitive Hernia Device Strength Evaluation*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2016. [Arbeitsplan]. WP108484.
15. Sharkawy AA, Klitzman B, Truskey GA, Reichert WM. Engineering the tissue which encapsulates subcutaneous implants. II. Plasma-tissue exchange properties. *Journal of Biomedical Materials Research* 1998;40(4):586–597.
16. Holt DJ, Grainger DW. Host response to biomaterials. In: Hollinger JO, ed. *An Introduction to Biomaterials*. 2nd ed. Boca Raton, FL: CRC Press; 2012;6:91–118.
17. Berman A. *Evaluation of Plexus without film and ETHICON PHYSIOMESH® in a 180-day rabbit subcutaneous model*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2015 [Studienprotokoll]. 23385C.
18. Berman A. *A Rabbit Model for the Biomechanical and Histological Assessment of Suture Line Wound Healing*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2008. [Studienprotokoll]. 19785C.
19. W. L. Gore & Associates, Inc. *GORE® BIO-A® Tissue Reinforcement Design Control (DC) Matrix*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2023. [Design Control Matrix – DC Matrix]. MD180587. 7. Fassung.
20. Morales-Conde S, Flores M, Fernandez V, Morales-Mendez S. *Bioabsorbable vs polypropylene plug for the “Mesh and Plug” inguinal hernia repair*. Poster präsentiert beim 9. Jahrestreffen der American Hernia Society. 9.-12. Februar 2005, San Diego, CA.
21. LeBlanc KA, Heniford BT, Voeller GR. Innovations in ventral hernia repair. Materials and techniques to reduce MRSA and other infections. *Contemporary Surgery* 2006;62(4)Supplement:1–8.
22. Clinger LR. *PTFE Knit Microbial Placement*. Flagstaff, AZ: W. L. Gore & Associates, Inc; 2018. [Arbeitsplan]. WP110158.
23. Lake SP, Ray S, Zihni AM, Thompson DM Jr, Gluckstein J, Deeken CR. Pore size and pore shape – but not mesh density – alter the mechanical strength of tissue ingrowth and host tissue response to synthetic mesh materials in a porcine model of ventral hernia repair. *Journal of the Mechanical Behavior of Biomedical Materials* 2015;42:186–197.
24. Klinge U, Klosterhalfen B, Birkenhauer V, Junge K, Conze J, Schumpelick V. Impact of polymer pore size on the interface scar formation in a rat model. *Journal of Surgical Research* 2002;103(2):208–214.
25. Zhu LM, Schuster P, Klinge U. Mesh implants: an overview of crucial mesh parameters. *World Journal of Gastrointestinal Surgery* 2015;7(10):226–236.
26. Klinge U, Klosterhalfen B, Conze J, et al. Modified mesh for hernia repair that is adapted to the physiology of the abdominal wall. *European Journal of Surgery* 1998;164(12):951–960.
27. Socea B, Socea LI, Bratu OG, et al. Recurrence rates and mesh shrinkage after polypropylene vs. polyester mesh hernia repair in complicated hernias. *Materiale Plastice* 2018;55(1):79–81.
28. Klosterhalfen B, Junge K, Klinge U. The lightweight and large porous mesh concept for hernia repair. *Expert Reviews in Medical Devices* 2005;2(1):103–117.
29. Jenkins ED, Yip M, Melman L, Frisella mm, Matthews BD. Informed consent: cultural and religious issues associated with the use of allogeneic and xenogeneic mesh products. *Journal of the American College of Surgeons* 2010;210(4):402–410.
30. Lighter M, Roberts J. The role of biosynthetic mesh in abdominal wall hernia repair in the setting of obesity, recurrence and high risk patients. Vorgestellt auf dem 2019 AHS Annual Meeting; 11.–14. März 2019; Las Vegas, NV. *Hernia* 2019;23 (Supplement 1): S96–S97. P-1256.
31. Landry M, Grimesly L, Gorman B, et al. Clinical Quality Improvement (CQI): evaluating the real work use of a new hybrid hernia mesh. Vorgestellt beim 2018 International Hernia Congress Program, 12. März 2018, Miami Beach, FL.
32. Kudsi OY, Gokcal F. Lateral approach totally extraperitoneal (TEP) robotic retromuscular ventral hernia repair. *Hernia* 2021;25(1):211–222.
33. W. L. Gore & Associates. *Assessment of GORE® SYNECOR Biomaterial in Focused Patient Populations and Long-Term Application*. NLM Identifier: NCT05094089. Veröffentlicht am 26. Oktober 2021. Aktualisiert am 20. April 2022. Letzter Zugriff am 19. Mai 2022. <https://clinicaltrials.gov/ct2/show/NCT05094089>
34. Yang XF, Liu JL. Anatomy essentials for laparoscopic inguinal hernia repair. *Annals of Translational Medicine* 2016;4(19):372.
35. Claus, C., Furtado, M., Malcher, F. et al. Ten golden rules for a safe MIS inguinal hernia repair using a new anatomical concept as a guide. *Surgical Endoscopy* 2020;34:1458–1464.

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© 2024 W. L. Gore & Associates GmbH 4GP2037-DE01 SEPTEMBER 2024
OPTILENE ist Markenzeichen von Aesculap AG.
BD, BARD und 3DMAX sind Markenzeichen von C.R. Bard, Inc.
TELA BIO und OVITEX sind Markenzeichen von TELA Bio, Inc.