

- ✓ No-profile radial support
- ✓ All-ePTFE unibody graft
- ✓ Cut and suture anywhere
- ✓ Smooth luminal surface



Smooth move.

Count on the leader in vascular graft technology to introduce a new option that offers all the benefits of a ringed graft—without the rings. This is like nothing you've seen or handled before. It's smooth to the touch but resists kinking and compression thanks to the all-ePTFE radial support concealed *within* the ultra-smooth graft wall. You can cut and suture with ease anywhere on the graft.

The patented unibody construction results in the lowest-profile radially supported graft available. It conforms exceptionally well to the native vessel and offers easy tunneling plus the potential for less tissue trauma. Of course, it also comes with Gore's heritage of product safety and excellence.

The GORE-TEX® INTERING vascular graft. Inner strength, outer elegance.



V A S C U L A R G R A F T



GORE, GORE-TEX®, INTERING, and designs are trademarks of W. L. Gore & Associates.
© 2001 W. L. Gore & Associates, Inc., Flagstaff, AZ 86004
AE0966-EN1 November 2001

www.goremedical.com EU: (00800)6334 4673

US: 800-437-8181