

Material matters when it comes to comfort!

GORE MYCROMESH® PLUS Biomaterial offers many benefits for soft tissue repair:

- ▶ Rapid tissue fixation
- ▶ Strong, effective repair
- ▶ Soft and conformable
- ▶ Trimmable
- ▶ Antimicrobial technology






MYCROMESH® PLUS

BIOMATERIAL

Composed of expanded polytetrafluoroethylene (ePTFE), one of the most inert polymers available, GORE MYCROMESH® PLUS Biomaterial features:

- Macroporous holes allow for rapid tissue fixation and incorporation**
 Experimental studies indicate the material is firmly adherent to host tissue after seven days with extensive vascularization (data on file).
- High material strength and suture retention**
 Regardless of the size or shape of the trimmed GORE MYCROMESH® PLUS Biomaterial, uniform strength is maintained.
- Material softness and conformability**
 Similar to host fascia, GORE MYCROMESH® PLUS Biomaterial is supple and non abrasive in an area of high flexion.
- Ease of handling**
 GORE MYCROMESH® PLUS Biomaterial may be trimmed and tailored without fraying.
- Effective against gram negative and gram positive organisms**
 Only Gore offers a soft tissue repair prosthesis with a broad-spectrum antimicrobial. GORE MYCROMESH® PLUS Biomaterial contains two antimicrobial agents which act synergistically to inhibit bacterial colonization of, and resist initial biofilm formation on, the material for up to 14 days post-implantation.
- Free of nerve entrapment**
 Unlike polypropylene, there are no reported hernia cases in the literature citing nerve entrapment when GORE MYCROMESH® PLUS Biomaterial was used.

Call us to evaluate the GORE MYCROMESH® Biomaterials in your next case.

GORE MYCROMESH® PLUS Biomaterial			Also available without antimicrobial technology GORE MYCROMESH® Biomaterial		
Catalogue Number	Nominal Thickness	Nominal Width x Length	Catalogue Number	Nominal Thickness	Nominal Width x Length
1MYMP01	1 mm	5 cm x 10 cm	1MYM01	1 mm	5 cm x 10 cm
1MYMP02	1 mm	5 cm x 15 cm	1MYM02	1 mm	5 cm x 15 cm
1MYMP06	1 mm	2.5 cm x 10 cm	1MYM07	1 mm	10 cm x 15 cm*
1MYMP07	1 mm	10 cm x 15 cm*	1MYM08	 1 mm	5 cm x 10 cm
1MYMP09	 1 mm	6 cm x 12 cm	1MYM09	 1 mm	6 cm x 12 cm
1MYMP10	1 mm	15 cm x 19 cm*	1MYM10	1 mm	15 cm x 19 cm*
1MYMP15	1 mm	18 cm x 24 cm	1MYM15	1 mm	18 cm x 24 cm
1MYMP16	1 mm	20 cm x 30 cm	1MYM16	1 mm	20 cm x 30 cm
1MYMP17	1 mm	26 cm x 34 cm*	1MYM17	1 mm	26 cm x 34 cm*

* Oval shaped

*Oval shaped



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

800.437.8181 (US)

928.779.2771 (US)

00800.6334.4673 (EU)

For international contacts and additional product information, visit goremedical.com