

Superior
TIPS
Performance



▶ **DEVICE SPECIFICATIONS**



VIATORR®

TIPS ENDOPROSTHESIS

PERFORMANCE
through data

GORE® VIATORR® TIPS Endoprosthesis

Designed for the Treatment of Portal Hypertension in Both *De Novo* and Revision Procedures

Catalogue Number	Endoprosthesis Internal Diameter (mm)	Graft-Lined Length (cm)	Graft Unlined Length (cm)	Maximum Guidewire Diameter (in)	Hemostatic Introducer Sheath (Fr)	Maximum Dilatation Balloon Diameter (mm)
PTB084275	8	4	2	≤ 0.038	10	8
PTB085275	8	5	2	≤ 0.038	10	8
PTB086275	8	6	2	≤ 0.038	10	8
PTB087275	8	7	2	≤ 0.038	10	8
PTB088275	8	8	2	≤ 0.038	10	8
PTB104275	10	4	2	≤ 0.038	10	10
PTB105275	10	5	2	≤ 0.038	10	10
PTB106275	10	6	2	≤ 0.038	10	10
PTB107275	10	7	2	≤ 0.038	10	10
PTB108275	10	8	2	≤ 0.038	10	10
PTB124275	12	4	2	≤ 0.038	10	12
PTB126275	12	6	2	≤ 0.038	10	12
PTB128275	12	8	2	≤ 0.038	10	12



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

 Consult Instructions for Use

Indications for Use: The GORE® VIATORR® TIPS Endoprosthesis is indicated for use in the *de novo* and revision treatment of portal hypertension and its complications such as variceal bleeding, gastropathy, refractory ascites, and / or hepatic hydrothorax. Refer to *Instructions for Use* at goremedical.com for a complete description of all contraindications, warnings, precautions and adverse events. Rx only

Products listed may not be available in all markets.

GORE®, PERFORMANCE THROUGH DATA, VIATORR®, and designs are trademarks of W. L. Gore & Associates. © 2004 – 2010 W. L. Gore & Associates, Inc. AH1501-EN2 AUGUST 2010