


CLINICAL REGISTRY SUMMARY

W.L. Gore & Associates conducted a Clinical Registry Program involving several centers in the US and one European center to support the safety of the GORE® ePTFE membranes in anterior spinal surgery.



GORE® ePTFE membranes were placed between vessels and various types of anterior spinal hardware to provide a permanent plane of dissection around vasculature to facilitate revisions.



- 478 documented implants
- ZERO reported device related complications

¹Anterior Lumbar Interbody Fusion
²Artificial Disc Replacement
³Anterior Cervical Discectomy and Fusion
⁴XLIF, Nucleus Replacement, Hybrid ADR/ALIF

Surgeon comments related to revisions reported at the time of data collection support the benefits of implanting GORE® PRECLUDE® Vessel Guard during the index procedure.

Contributors:

Gary R. Bench, MD, San Diego General Surgery, San Diego, CA
Rudolph Bertagnoli, MD, Spine Center Straubing, Bogen, Germany
Salvador Brau, MD, FACS, Spine Access Surgery Associates, Los Angeles, CA
Earl D. Cottrell, MD, General Vascular Specialists, Las Vegas, NV
William Lu, MD, Orlando Neurosurgery, Winter Park, FL
James G. Melton, DO, Oklahoma Cardiovascular Associates, Oklahoma City, OK
Thomas Terramani, MD, Vascular Associates of San Diego, San Diego, CA

W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com



V E S S E L G U A R D

Products listed may not be available in all markets.
GORE®, PERFORMANCE THROUGH EXPERIENCE, PRECLUDE®, and designs are trademarks of W. L. Gore & Associates.
© 2007, 2008, 2010 W. L. Gore & Associates, Inc. AL2903-EN5 MARCH 2010

PERFORMANCE
through experience