



GORE Tri-Lobe Balloon Catheter Inflation Volume Guide



TRI-LOBE

BALLOON CATHETER



INTENDED BALLOON DIAMETER (mm)	RECOMMENDED INFLATION VOLUME (cc)		
	Inner Vessel Diameter (mm)	BCM1634	BCL2645
16		5	
17		5	
18		5	
19		6	
20		6	
21		6	
22		7	
23		8	
24		8	
25		9	
26		10	11
27		11	12
28		11	12
29		12	13
30		13	14
31		14	15
32		15	16
33		16	18
34		17	19
35			19
36			20
37			21
38			22
39			24
40			25
41			27
42			28

Products listed may not be available in all markets.
 GORE and designs are trademarks of W. L. Gore & Associates.
 © 2010 W. L. Gore & Associates, Inc. AP0128-EN1 MARCH 2010