



## CONFIGURATIONS

Catalogue number	Description	Treatment range
ASD27E	27 mm GORE® CARDIOFORM ASD Occluder	8–15 mm
ASD32E	32 mm GORE® CARDIOFORM ASD Occluder	13–20 mm
ASD37E	37 mm GORE® CARDIOFORM ASD Occluder	18–25 mm
ASD44E	44 mm GORE® CARDIOFORM ASD Occluder	23–30 mm
ASD48E	48 mm GORE® CARDIOFORM ASD Occluder	28–35 mm

All device configurations have an 80 cm working length and 10–14 Fr working catheter depending on occluder size.

Refer to *Instructions for Use* at [goremedical.com](http://goremedical.com) for a complete description of all contraindications, warnings, precautions, and adverse events.

Products listed may not be available in all markets.

GORE, *Together, improving life*, CARDIOFORM and designs are trademarks of W. L. Gore & Associates, Inc.  
© 2019 W. L. Gore & Associates GmbH JULY 2019

---

**W. L. Gore & Associates, Inc.**  
Flagstaff, AZ 86004  
+65.67332882 (Asia Pacific) 800.437.8181 (United States)  
00800.6334.4673 (Europe) 928.779.2771 (United States) [goremedical.com](http://goremedical.com)

