


ORDERING INFORMATION

Needle Coding System



THX-26

Number describes the length
of the needle in mm

First letter describes
needle type
T = Taper Point
P = Piercing Point
R = Reverse Cutting
C = Cutting

Second letter describes
curvature of the needle
L = 1/2 curve (ski)
T = 3/8 circle
H = 1/2 circle
F = 5/8 circle
S = Straight

Third letter, if present,
denotes the following:
X = Heavier wire diameter, e.g. THX-36
c = 'c' series, premium needle

EXAMPLE

TTC-13 = Taper Point, 3/8 circle,
premium needle,
13 mm in length

PH-17 = Piercing Point, 1/2 circle,
17 mm in length

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. ℞ only

Products listed may not be available in all markets.

GORE, GORE-TEX, *Together, improving life* and designs are trademarks of W. L. Gore & Associates.

© 1996, 1998, 2000, 2001, 2005, 2007, 2008, 2012, 2014, 2023, 2025 W. L. Gore & Associates GmbH 25PL2018-EN01 JANUARY 2025

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 6733 2882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771



				Thread Size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	CV-0
				Approx USP Size:	7-0	6-0	5-0	4-0	3-0	2-0	0	1
	Single-armed Needles	Double-armed Needles	Needle Code	Thread Length	Catalogue Number							
○ Taper Point $\frac{3}{8}$ Circle			TTC-7	24" – 61 cm	8K10E							
			TTC-9	18" – 46 cm	8J02E	7J02E	6J02E					
				24" – 61 cm	8K02E	7K02E	6K02E					
				30" – 76 cm	8M02E	7M02E	6M02E					
			TTC-12	30" – 76 cm		7M12E	6M10E					
			TTC-13	24" – 61 cm		7K06E		5K02E				
				30" – 76 cm	8M04E	7M06E	6M04E	5M02E				
				36" – 91 cm			6N04E	5N02E				
			TT-17	36" – 91 cm				5N06E	4N04E	3N02E		
○ Taper Point $\frac{1}{2}$ Circle			THC-9	30" – 76 cm			6M24E					
			THC-13	24" – 61 cm				5K12E				
				30" – 76 cm		7M10E	6M12E	5M12E				
				36" – 91 cm				5N12E				
			TH-18	36" – 91 cm			6N16E	5N04E	4N02E	3N04E		
				48" – 122 cm					4U02E			
			TH-26	36" – 91 cm				5N20E	4N06E	3N08E	2N02E	
				36" – 91 cm							2N01E	
				48" – 122 cm					4U06E	3U08E		
			TH-35	36" – 91 cm					4N12E	3N06E		
				48" – 122 cm						3U06E		
			TH-50	48" – 122 cm								0U01E
			THX-22	36" – 91 cm							2N03E	
			THX-25	36" – 91 cm							2N13E	0N13E
			THX-26	36" – 91 cm							2N06E	
				36" – 91 cm							2N05E	0N05E
				48" – 122 cm							0U06E	
			THX-36	36" – 91 cm							2N08E	0N08E
				36" – 91 cm							2N07E	0N07E
○ Taper Point $\frac{1}{4}$ Circle			TFX-25	36" – 91 cm							2N19E	0N09E
○ Taper Point $\frac{1}{2}$ Curve (ski) Smooth, even point for easily penetrated tissues, tapered body			TLX-25	48" – 122 cm								0U03E
○ Taper Point Straight			TS-50	36" – 91 cm							2N23E	
▲ Piercing Point $\frac{3}{8}$ Circle Multifaceted tip for tough tissue, tapered body			PT-9	24" – 61 cm	8K08E	7K04E	6K06E					
				30" – 76 cm	8M08E	7M04E	6M06E					
			PT-13	30" – 76 cm	8M06E	7M08E	6M08E					
				36" – 91 cm				5N08E				
			PT-17	36" – 91 cm				5N10E	4N10E	3N14E		
▲ Piercing Point $\frac{1}{2}$ Circle			PH-13	30" – 76 cm			6M14E	5M16E				
			PH-17	36" – 91 cm				5N14E	4N08E	3N12E		
			PH-24	36" – 91 cm					4N16E	3N10E	2N10E	
				48" – 122 cm						3U10E		
			PH-35	48" – 122 cm						3U20E		
▲ Straight Cutting Cutting edges for tough tissue, triangular body			CS-65	30" – 76 cm							2M25E	

Catalogue number in bold type indicates 1:1 needle to thread ratio to minimize suture hole bleeding.

Telephone Orders:
800.437.8181, 928.779.2771 (US)
00800.6334.4673 (EU)

