






OPEN MORE POSSIBILITIES: PREDICTABLE PERFORMANCE WITH ENHANCEMENTS IN ACCESS AND DELIVERY

Large diameter configurations with lower delivery profile

- Expanded accessibility through smaller sheaths.
- Smaller sheath sizes result in a lower risk of vascular access complications in select patient populations.¹
- Enhanced visualization under fluoroscopy
 - Four gold radiopaque markers on the distal and proximal ends for more confident delivery in complex cases



Endoprosthesis labeled diameter (mm)*	Recommended vessel diameter (mm)†	Legacy profile (Fr)	New lower profile (Fr)‡
9	7.6–8.5	9	8 
10	8.6–9.5	11	8 
11	9.6–10.5	11	10 
13	10.6–12.0	12	10† 

* Labeled device diameters and lengths are nominal.


† Recommended endoprosthesis compression within the vessel is approximately 5–20%.

‡ The 13 mm diameter device is not compatible with the 10 Fr COOK® FLEXOR® CHECK-FLO® Sheath.

§ Reduced profile compared to prior size. Not to scale.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

Continued innovation for durable outcomes and unmatched versatility



* GORE® VIABAHN® Endoprosthesis. W. L. Gore & Associates website. Accessed July 29, 2020. <https://www.goremedical.com/VIABAHN/references>.

Reference

1. Applegate RJ, Sacrinty MT, Kutcher MA, et al. Trends in vascular complications after diagnostic cardiac catheterization and percutaneous coronary intervention via the femoral artery, 1998 to 2007. *JACC: Cardiovascular Interventions* 2008;1(3):317–326.

 Consult Instructions for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

Products listed may not be available in all markets.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

GORE, *Together, improving life*, HEMOBAHN, PROPATEN, VIABAHN and designs are trademarks of W. L. Gore & Associates.

© 2021 W. L. Gore & Associates GmbH 2157883-EN1 APRIL 2021

W. L. Gore & Associates, Inc.

goremedical.com

Asia Pacific +65 6733 2882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673

United States Flagstaff, AZ 86003 800 437 8181 928 779 2771

