

GORE® VIABAHN® VBX
Balloon Expandable Endoprosthesis



NOW
6Fr
compatible
**GREATER
VERSATILITY**

CONFIGURATIONS

Together, improving life



Catalogue numbers

GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis (VBX Stent Graft)

VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXB051502E	5	8	15	135	.035	6
BXB051902E	5	8	19	135	.035	6
BXB052902E	5	8	29	135	.035	6
BXB053902E	5	8	39	135	.035	6
BXB055902E	5	8	59	135	.035	6
BXB057902E	5	8	79	135	.035	6
BXB061502E	6	8	15	135	.035	6
BXB061902E	6	8	19	135	.035	6
BXB062902E	6	8	29	135	.035	6
BXB063902E	6	8	39	135	.035	6
BXB065902E	6	8	59	135	.035	6
BXB067902E	6	8	79	135	.035	6
BXB071502E	7	11	15	135	.035	6
BXB071902E	7	11	19	135	.035	6
BXB072902E	7	11	29	135	.035	6
BXB073902E	7	11	39	135	.035	6
BXB075902E	7	11	59	135	.035	6
BXB077902E	7	11	79	135	.035	6
BXB082902E	8	11	29	135	.035	7
BXB083902E	8	11	39	135	.035	7
BXB085902E	8	11	59	135	.035	7
BXB087902E	8	11	79	135	.035	7
BXB092902E	9	13	29	135	.035	7
BXB093902E	9	13	39	135	.035	7
BXB095902E	9	13	59	135	.035	7

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXB097902E	9	13	79	135	.035	7
BXA102902E†	10	13	29	135	.035	8
BXA103902E†	10	13	39	135	.035	8
BXA105902E†	10	13	59	135	.035	8
BXA107902E†	10	13	79	135	.035	8
BXA112902E†	11	16	29	135	.035	8
BXA113902E†	11	16	39	135	.035	8
BXA115902E†	11	16	59	135	.035	8
BXA117902E†	11	16	79	135	.035	8
BXBL082902E	8	16	29	135	.035	7
BXBL083902E	8	16	39	135	.035	7
BXBL085902E	8	16	59	135	.035	8
BXBL087902E	8	16	79	135	.035	8

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

† The following larger diameter devices are unchanged and will continue to use the same catalogue number.

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx only
Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2024 W. L. Gore & Associates, Inc. 24PL2053-EN01 JUNE 2024

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

