



PREDICTABLE PERFORMANCE

With lower profile, radiopaque markers and the CBAS Heparin Surface

Lower delivery profile on large diameter configurations

- Expanded accessibility through smaller sheaths for 9–13 mm diameter devices
- Smaller sheath sizes result in a lower risk of vascular access complications in select patient populations¹

Endoprosthesis labeled diameter (mm) [†]	Recommended vessel diameter (mm) [‡]	Legacy profile (Fr)	New lower profile (Fr)
9	7.6–8.5	○ 9	○ 8
10	8.6–9.5	○ 11	○ 8
11	9.6–10.5	○ 11	○ 10
13	10.6–12.0	○ 12	○ 10 [§]

* As used by Gore, Heparin Bioactive Surface refers to Gore's proprietary CBAS Heparin Surface.

[†] Labeled device diameters and lengths are nominal.

[‡] Recommended endoprosthesis compression within the vessel is approximately 5–20%.

[§] The 13 mm diameter device is not compatible with the 10 Fr COOK® FLEXOR® CHECK-FLO® Introducers.


COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface

Continued innovation for durable outcomes and unmatched versatility


Enhanced visualization under fluoroscopy

- Four gold radiopaque markers on the distal and proximal ends of all device configurations for more confident delivery



Lasting thromboresistance, proven technology²


- CBAS Heparin Surface is the proven heparin bonding technology for lasting thromboresistance
- Proven heparin availability and bioactivity for an unmatched and persistent ability to take up antithrombin
- See the technology at [goremedical.com/la/cbas](https://www.goremedical.com/la/cbas)



References

- Applegate RJ, Sacrinty MT, Kutcher MA, *et al.* Trends in vascular complications after diagnostic cardiac catheterization and percutaneous coronary intervention via the femoral artery, 1998 to 2007. *JACC: Cardiovascular Interventions* 2008;1(3):317-326.
- CBAS Heparin Surface. W. L. Gore & Associates website. Accessed October 28, 2020. <https://www.goremedical.com/cbas/references>

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. 

Products listed may not be available in all markets.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.
GORE, *Together, improving life*, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2021 W. L. Gore & Associates, Inc. 2139516-EN FEBRUARY 2021

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

