

ORDERING INFORMATION





Needle-to-thread ratio

Thread size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	CV-0
Approximate USP size:	7-0	6-0	5-0	4-0	3-0	2-0	0	1
	TTc-7 (1.6:1)	TTc-9 (1.4:1)	TTc-9 (1:1)	TTc-13 (1:1)	TT-17 (1.2:1)	TT-17 (1:1)	TFX-25 (1.6:1)	TFX-25 (1.6:1)
	TTc-9 (1.9:1)	TTc-12 (1.4:1)	TTc-12 (1:1)	TT-17 (1.2:1)	TH-18 (1.2:1)	TH-18 (1:1)	TH-26 (1:1)	TLX-25 (1.6:1)
	TTc-13 (2.2:1)	TTc-13 (1.6:1)	TTc-13 (1.4:1)	THc-13 (1:1)	TH-26 (1.6:1)	TH-26 (1.2:1)	THX-22 (1.9:1)	TH-50 (1.6:1)
	PT-9 (1.9:1)	THc-13 (1.6:1)	THc-9 (1:1)	TH-18 (1.2:1)	TH-35 (1.6:1)	TH-35 (1.2:1)	THX-25 (1.9:1)	THX-25 (1.6:1)
	PT-13 (2.2:1)	PT-9 (1.4:1)	THc-13 (1.4:1)	TH-26 (1.7:1)	PT-17 (1.2:1)	PT-17 (1:1)	THX-26 (1.6:1)	THX-26 (1.3:1)
		PT-13 (1.6:1)	TH-18 (1.6:1)	PT-13 (1:1)	PH-17 (1.2:1)	PH-17 (1:1)	THX-36 (1.6:1)	THX-36 (1.3:1)
			PT-9 (1:1)	PT-17 (1.2:1)	PH-24 (1.6:1)	PH-24 (1.2:1)	PH-24 (1:1)	
			PT-13 (1.4:1)	PH-13 (1:1)		PH-30 (1.3:1)	CS-65 (1.2:1)	
			PH-13 (1.4:1)	PH-17 (1.2:1)		PH-35 (1.5:1)	TS-50 (1.2:1)	

Indicates products with 1:1 needle-to-thread ratio

Needle coding system



THX-26

Number describes the length of the needle in mm







First letter describes needle type
T = Taper point
P = Piercing point
R = Reverse cutting
C = Cutting

Second letter describes curvature of the needle
L = 1/2 curve (ski)
T = 3/8 circle
H = 1/2 circle
F = 5/8 circle
S = Straight

Third letter, if present, denotes the following:
X = Heavier wire diameter, e.g. THX-36 c = 'c' series, premium needle

Example
TTc-13 = Taper point, 3/8 circle, premium needle, 13 mm in length
PH-17 = Piercing point, 1/2 circle, 17 mm in length



		Thread size:															
		Approx USP size:															
		CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	CV-0								
		7-0	6-0	5-0	4-0	3-0	2-0	0	1								
		Single-armed needles	Double-armed needles	Needle code	Thread length	Catalogue number								Ethicon reference	Medtronic reference		
⦿	Taper point ⅜ circle			TTc-7	24" – 61 cm	8K10								BV130-5	CV-351, MV-135-5		
				TTc-9	18" – 46 cm	8J02	7J02	6J02							BV-1, BV175-8	CV-1, CVF-1, MV-175-8	
			24" – 61 cm		8K02	7K02	6K02										
			30" – 76 cm		8M02	7M02	6M02										
					TTc-12	30" – 76 cm		7M12	6M10						BV	CV-11, CVF	
					TTc-13	24" – 61 cm		7K06		5K02						C-1	CV-11, CVF-11
		30" – 76 cm	8M04	7M06		6M04	5M02										
		36" – 91 cm				6N04	5N02										
			TT-17	36" – 91 cm				5N06	4N04	3N02			BB	CV-15, CVF-15			
⦿	Taper point ½ circle			THc-9	30" – 76 cm			6M24						N/A	CV-20		
				THc-13	18" – 46 cm		7J10								RB-2	CV-22, CVF-22	
			24" – 61 cm					5K12									
			30" – 76 cm			7M10	6M12	5M12									
			36" – 91 cm					5N12									
					TH-18	36" – 91 cm			6N16	5N04	4N02	3N04			RB-1	CV-23, CVF-23	
			48" – 122 cm						4U02								
				TH-26	36" – 91 cm				5N20	4N06	3N08	2N02			SH	V-20	
			36" – 91 cm							2N01							
			48" – 122 cm					4U06	3U08								
					TH-35	36" – 91 cm					4N12	3N06			MH	V-26	
			48" – 122 cm							3U06							
				TH-50	48" – 122 cm									OU01	CTX	GS-25	
				THX-22	36" – 91 cm								2N03		M07	GS-23	
				THX-25	36" – 91 cm								2N13	0N13	M06	GS-22	
		THX-26	36" – 91 cm								2N06		CT-2	GS-22			
			36" – 91 cm							2N05	0N05						
			48" – 122 cm							OU06							
		THX-36	36" – 91 cm								2N08	0N08	CT-1	GS-21			
			36" – 91 cm							2N07	0N07						
⦿	Taper point ⅜ circle			TFX-25	36" – 91 cm							2N19	0N09	UR-6	GU-46		
⦿	Taper point ½ curve (ski) Smooth, even point for easily penetrated tissues, tapered body			TLX-25	48" – 122 cm								OU03	EN-5	SK		
⦿	Taper point straight			TS-50	36" – 91 cm							2N23		ST-1	ST-1		
⬆	Piercing point ⅜ circle Multifaceted tip for tough tissue, tapered body			PT-9	24" – 61 cm	8K08	7K04	6K06						CC	KV-1		
					30" – 76 cm	8M08	7M04	6M06									
				PT-13	30" – 76 cm	8M06	7M08	6M08						CC-1	KV-11		
				PT-17	36" – 91 cm				5N10	4N10	3N14			V-4	KV-15		
⬆	Piercing point ½ circle			PH-13	30" – 76 cm			6M14	5M16					N/A	N/A		
				PH-17	36" – 91 cm				5N14	4N08	3N12			V-5	KV-5		
Multi-needle		PH-24	36" – 91 cm					4N16	3N10	2N10		V-7	KV-7				
			48" – 122 cm						3U10								
			PH-35	48" – 122 cm						3U20			V-34	KV-34			
⬆	Straight cutting Cutting edges for tough tissue, triangular body			CS-65	30" – 76 cm							2M25		KS	N/A		

Catalogue number in bold type indicates 1:1 needle-to-thread ratio to minimize suture hole bleeding.

To order a one-dozen box, suffix the catalogue number with the letter A.

Telephone orders:
 U.S. toll free: +1 800 528 8763
 U.S. direct: +1 928 864 2927
 Europe: +00800 6334 4673

 Consult Instructions
for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. **R_X Only**

Products listed may not be available in all markets.

GORE, GORE-TEX, *Together, improving life* and designs are trademarks of W. L. Gore & Associates.
 © 1996, 1998, 2000, 2001, 2005, 2007, 2008, 2012, 2014, 2023, 2025 W. L. Gore & Associates, Inc. 24SS1007-EN01 APRIL 2025

W. L. Gore & Associates, Inc.
 goremedical.com
Asia Pacific +65 6733 2882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

