


GORE® VIABAHN® VBX
Balloon Expandable Endoprosthesis



NOW
6Fr
compatible
**GREATER
VERSATILITY**

CONFIGURATIONS

Together, improving life



Catalogue numbers

GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis (VBX Stent Graft)

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/ nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXA051502B	BXB051502B	5	8	15	135	.035	6
BXA051902B	BXB051902B	5	8	19	135	.035	6
BXA052902B	BXB052902B	5	8	29	135	.035	6
BXA053902B	BXB053902B	5	8	39	135	.035	6
BXA055902B	BXB055902B	5	8	59	135	.035	6
BXA057902B	BXB057902B	5	8	79	135	.035	6
BXA061502B	BXB061502B	6	8	15	135	.035	6
BXA061902B	BXB061902B	6	8	19	135	.035	6
BXA062902B	BXB062902B	6	8	29	135	.035	6
BXA063902B	BXB063902B	6	8	39	135	.035	6
BXA065902B	BXB065902B	6	8	59	135	.035	6
BXA067902B	BXB067902B	6	8	79	135	.035	6
BXA071502B	BXB071502B	7	11	15	135	.035	6
BXA071902B	BXB071902B	7	11	19	135	.035	6
BXA072902B	BXB072902B	7	11	29	135	.035	6
BXA073902B	BXB073902B	7	11	39	135	.035	6
BXA075902B	BXB075902B	7	11	59	135	.035	6
BXA077902B	BXB077902B	7	11	79	135	.035	6
BXA082902B	BXB082902B	8	11	29	135	.035	7
BXA083902B	BXB083902B	8	11	39	135	.035	7
BXA085902B	BXB085902B	8	11	59	135	.035	7
BXA087902B	BXB087902B	8	11	79	135	.035	7

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXA092902B	BXB092902B	9	13	29	135	.035	7
BXA093902B	BXB093902B	9	13	39	135	.035	7
BXA095902B	BXB095902B	9	13	59	135	.035	7
BXA097902B	BXB097902B	9	13	79	135	.035	7
BXA102902B	Unchanged	10	13	29	135	.035	8
BXA103902B	Unchanged	10	13	39	135	.035	8
BXA105902B	Unchanged	10	13	59	135	.035	8
BXA107902B	Unchanged	10	13	79	135	.035	8
BXA112902B	Unchanged	11	16	29	135	.035	8
BXA113902B	Unchanged	11	16	39	135	.035	8
BXA115902B	Unchanged	11	16	59	135	.035	8
BXA117902B	Unchanged	11	16	79	135	.035	8
BXAL082902B	BXBL082902B	8	16	29	135	.035	7
BXAL083902B	BXBL083902B	8	16	39	135	.035	7
BXAL085902B	BXBL085902B	8	16	59	135	.035	8
BXAL087902B	BXBL087902B	8	16	79	135	.035	8

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).



Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. R_x Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2023–2025 W. L. Gore & Associates, Inc. 24PL1191-EN01 JANUARY 2025

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 **Australia/New Zealand** 1800 680 424 **Europe** 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

