


GORE® VIABAHN® VBX
Balloon Expandable Endoprosthesis



NOW
6Fr
compatible
**GREATER
VERSATILITY**

CONFIGURATIONS

Together, improving life



Catalogue numbers

GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis (VBX Stent Graft)

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/ nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXA051501A	BXB051501A	5	8	15	80	.035	6
BXA051901A	BXB051901A	5	8	19	80	.035	6
BXA052901A	BXB052901A	5	8	29	80	.035	6
BXA053901A	BXB053901A	5	8	39	80	.035	6
BXA055901A	BXB055901A	5	8	59	80	.035	6
BXA057901A	BXB057901A	5	8	79	80	.035	6
BXA061501A	BXB061501A	6	8	15	80	.035	6
BXA061901A	BXB061901A	6	8	19	80	.035	6
BXA062901A	BXB062901A	6	8	29	80	.035	6
BXA063901A	BXB063901A	6	8	39	80	.035	6
BXA065901A	BXB065901A	6	8	59	80	.035	6
BXA067901A	BXB067901A	6	8	79	80	.035	6
BXA071501A	BXB071501A	7	11	15	80	.035	6
BXA071901A	BXB071901A	7	11	19	80	.035	6
BXA072901A	BXB072901A	7	11	29	80	.035	6
BXA073901A	BXB073901A	7	11	39	80	.035	6
BXA075901A	BXB075901A	7	11	59	80	.035	6
BXA077901A	BXB077901A	7	11	79	80	.035	6
BXA082901A	BXB082901A	8	11	29	80	.035	7
BXA083901A	BXB083901A	8	11	39	80	.035	7
BXA085901A	BXB085901A	8	11	59	80	.035	7
BXA087901A	BXB087901A	8	11	79	80	.035	7

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/ nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXA092901A	BXB092901A	9	13	29	80	.035	7
BXA093901A	BXB093901A	9	13	39	80	.035	7
BXA095901A	BXB095901A	9	13	59	80	.035	7
BXA097901A	BXB097901A	9	13	79	80	.035	7
BXA102901A	Unchanged	10	13	29	80	.035	8
BXA103901A	Unchanged	10	13	39	80	.035	8
BXA105901A	Unchanged	10	13	59	80	.035	8
BXA107901A	Unchanged	10	13	79	80	.035	8
BXA112901A	Unchanged	11	16 [†]	29	80	.035	8
BXA113901A	Unchanged	11	16 [†]	39	80	.035	8
BXA115901A	Unchanged	11	16 [†]	59	80	.035	8
BXA117901A	Unchanged	11	16 [†]	79	80	.035	8
BXAL082901A	BXBL082901A	8	16 [†]	29	80	.035	7
BXAL083901A	BXBL083901A	8	16 [†]	39	80	.035	7
BXAL085901A	BXBL085901A	8	16 [†]	59	80	.035	8
BXAL0857901A	BXBL0857901A	8	16 [†]	79	80	.035	8
BXA051502A	BXB051502A	5	8	15	135	.035	6
BXA051902A	BXB051902A	5	8	19	135	.035	6
BXA052902A	BXB052902A	5	8	29	135	.035	6
BXA053902A	BXB053902A	5	8	39	135	.035	6
BXA055902A	BXB055902A	5	8	59	135	.035	6
BXA057902A	BXB057902A	5	8	79	135	.035	6
BXA061502A	BXB061502A	6	8	15	135	.035	6
BXA061902A	BXB061902A	6	8	19	135	.035	6
BXA062902A	BXB062902A	6	8	29	135	.035	6
BXA063902A	BXB063902A	6	8	39	135	.035	6
BXA065902A	BXB065902A	6	8	59	135	.035	6
BXA067902A	BXB067902A	6	8	79	135	.035	6

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

† Technical limit of the device as determined by in-vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the VBX FLEX Clinical Study and is outside of the approved indication – see *Instructions for Use*.

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm)*	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Introducer sheath size (Fr)
BXA071502A	BXB071502A	7	11	15	135	.035	6
BXA071902A	BXB071902A	7	11	19	135	.035	6
BXA072902A	BXB072902A	7	11	29	135	.035	6
BXA073902A	BXB073902A	7	11	39	135	.035	6
BXA075902A	BXB075902A	7	11	59	135	.035	6
BXA077902A	BXB077902A	7	11	79	135	.035	6
BXA082902A	BXB082902A	8	11	29	135	.035	7
BXA083902A	BXB083902A	8	11	39	135	.035	7
BXA085902A	BXB085902A	8	11	59	135	.035	7
BXA087902A	BXB087902A	8	11	79	135	.035	7
BXA092902A	BXB092902A	9	13	29	135	.035	7
BXA093902A	BXB093902A	9	13	39	135	.035	7
BXA095902A	BXB095902A	9	13	59	135	.035	7
BXA097902A	BXB097902A	9	13	79	135	.035	7
BXA102902A	Unchanged	10	13	29	135	.035	8
BXA103902A	Unchanged	10	13	39	135	.035	8
BXA105902A	Unchanged	10	13	59	135	.035	8
BXA107902A	Unchanged	10	13	79	135	.035	8
BXA112902A	Unchanged	11	16 [†]	29	135	.035	8
BXA113902A	Unchanged	11	16 [†]	39	135	.035	8
BXA115902A	Unchanged	11	16 [†]	59	135	.035	8
BXA117902A	Unchanged	11	16 [†]	79	135	.035	8
BXAL082902A	BXBL082902A	8	16 [†]	29	135	.035	7
BXAL083902A	BXBL083902A	8	16 [†]	39	135	.035	7
BXAL085902A	BXBL085902A	8	16 [†]	59	135	.035	8
BXAL087902A	BXBL087902A	8	16 [†]	79	135	.035	8

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (Secondary balloon not included).

† Technical limit of the device as determined by in-vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the VBX FLEX Clinical Study and is outside of the approved indication – see *Instructions for Use*.

 Consult Instructions
for Use
eifu.goremedical.com

INDICATIONS FOR USE IN THE U.S.: The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis is indicated for the treatment of de novo or restenotic lesions found in iliac arteries with reference vessel diameters ranging from 5 mm–13 mm and lesion lengths up to 110 mm, including lesions at the aortic bifurcation. **CONTRAINDICATIONS:** Do not use the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis in patients with known hypersensitivity to heparin, including those patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II. Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. R_x Only

Products listed may not be available in all markets.

GORE, *Together, improving life*, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2023 W. L. Gore & Associates, Inc. 231095892-EN NOVEMBER 2023

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

