


GORE-TEX® SUTURE: THE PERFECT CLOSE TO YOUR VASCULAR PROCEDURES

			Thread size:				CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	
			Approx USP size:				7-0	6-0	5-0	4-0	3-0	2-0	0	
	Double-armed needles	Needle code	JOHNSON & JOHNSON® ETHICON® needle reference	MEDTRONIC SYNETURE needle reference	Thread length	Catalogue number								
Taper point 3/8 circle 		TTC-7	BV130-5	CV-351, MV-135-5	24" – 61 cm	8K10								
		TTC-9	BV-1, BV175-8	CV-1, CVF-1, MV-175-8	18" – 46 cm	8J02	7J02	6J02						
					24" – 61 cm	8K02	7K02	6K02						
		TTC-12	BV	CV-11, CVF	30" – 76 cm	8M02	7M02	6M02						
					24" – 61 cm		7M12	6K10						
		TTC-13	C-1	CV-11, CVF-11	30" – 76 cm		7J06	6J04	6M10					
24" – 61 cm						7K06	6K04	5K02						
Taper point 1/2 circle 		TT-17	BB	CV-15, CVF-15	36" – 91 cm				6M04	6M04	5M02			
								6N04	5N02					
		THC-9	N/A	CV-20	30" – 76 cm				5N06					
									5N06	4N04	3N02			
		THC-13	RB-2	CV-22, CVF-22	18" – 46 cm		7J10	6J12						
					24" – 61 cm			6K12	5K12					
	TH-18	RB-1	CV-23, CVF-23	30" – 76 cm		7M10	6M12	5M12						
				36" – 91 cm				5N12						
Piercing point 3/8 circle Multifaceted tip for tough tissue, tapered body 		TH-22 (with pledgets)	SH-1	CV-25	36" – 91 cm				6N16	5N04	4N02	3N04		
					48" – 122 cm					4U02				
		TH-26	SH	V-20	36" – 91 cm					3201*	3200*			
					48" – 122 cm									
		TH-35	MH	V-26	36" – 91 cm					5N20	4N06	3N08	2N02	
					48" – 122 cm					4U06	3U08	2U02		
Piercing point 1/2 circle Multifaceted tip for tough tissue, tapered body 		PT-9	CC	KV-1	18" – 46 cm		7J04							
					24" – 61 cm	8K08	7K04	6K06						
		PT-13	CC-1	KV-11	30" – 76 cm	8M08	7M04	6M06						
					24" – 61 cm		7K08	6K08	5K08					
		PT-17	V-4	KV-15	30" – 76 cm	8M06	7M08	6M08						
					36" – 91 cm				5N08					
Piercing point 3/8 circle Multifaceted tip for tough tissue, tapered body 		PT-17	V-4	KV-15	36" – 91 cm				5N10	4N10	3N14			

Catalogue numbers in bold type indicate 1:1 needle to thread ratio to minimize suture hole bleeding.

JOHNSON & JOHNSON and ETHICON are trademarks of Johnson & Johnson Company.
MEDTRONIC and SYNETURE are trademarks of Medtronic, Inc.

		Thread size:		CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2
		Approx USP size:		7-0	6-0	5-0	4-0	3-0	2-0	0
	Double-armed needles	Needle code	JOHNSON & JOHNSON® ETHICON® needle reference	MEDTRONIC SYNETURE needle reference	Thread length	Catalogue number				
Piercing point 1/2 circle Ⓐ		PH-13	N/A	N/A	30" – 76 cm	6M14	5M16			
		PH-17	V-5	KV-5	36" – 91 cm		5N14	4N08	3N12	
		PH-24	V-7	KV-7	36" – 91 cm			4N16	3N10	2N10
					48" – 122 cm			4U16	3U10	2U10
		PH-30	N/A	N/A	42" – 107 cm				3P24	
		PH-35	V-34	KV-34	48" – 122 cm				3U20	

Catalogue numbers in bold type indicate 1:1 needle to thread ratio to minimize suture hole bleeding.

To order a three-dozen box, suffix the catalogue number with the letter B. To order a one-dozen box, suffix the catalogue number with the letter A. Some catalogue numbers available in A (one-dozen) boxes only.

Commonly requested GORE-TEX® Sutures

Thread size	Thread length (cm)	Taper point needles	Catalogue number	Piercing point needles	Catalogue number
Proximal aorto-bifemoral					
CV-3	91	TH-26	3N08	PH-24	3N10
CV-4	91	TH-26	4N06	PH-24	4N16
	91	TH-18	4N02		
Distal aorto-bifemoral					
	76	THC-13	5M12		
CV-5	76	THC-13	5M02		
	91	TH-18	5N04		
CV-6	76	TTc-12	6M10	PT-13	6M08
Proximal femoro-popliteal					
CV-5	91	TTc-13	5N02	PT-13	5N08
CV-6	76	TTc-9	6M02		
CV-6	76	TTc-12	6M10		
CV-6	76	TTc-13	6M04	PT-13	6M08
Distal femoro-popliteal					
CV-6	76	TTc-9	6M02	PT-9	6M06
CV-7	61	TTc-9	7K02		
	76	TTc-9	7K02	PT-9	6M06
CV-8	61	TTc-9	8K02		

Thread size	Thread length (cm)	Taper point needles	Catalogue number	Piercing point needles	Catalogue number
Carotid endarterectomy					
CV-6	76	TTc-9	6M02		
	76	THC-13	6M04		
Axillo-femoral					
CV-5	91	TTc-13	5N02	PT-13	5N08
	76	TTc-9	6M02		
CV-6	76	TTc-12	6M10		
	76	TTc-13	6M04	PT-13	6M08
AV access graft					
	61/76	TTc-9	6K02/6M02		
CV-6	61/76	TTc-12	6K10/6M10		
	61/76	TTc-13	6K04/6M04		
CV-7	61/76	TTc-9	7K02/7M02		
AV fistulas					
CV-7	46	TTc-9	7J02		
CV-8	46	TTc-9	8J02		

Catalogue numbers in bold type indicate 1:1 needle to thread ratio to minimize suture hole bleeding.

* With ePTFE pledgets (3).

For a full listing of available GORE-TEX® Suture, visit goremedical.com/products/suture



Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. 

Products listed may not be available in all markets.

JOHNSON & JOHNSON and ETHICON are trademarks of Johnson & Johnson Company. MEDTRONIC and SYNETURE are trademarks of Medtronic, Inc.

GORE, GORE-TEX, *Together, improving life* and designs are trademarks of W. L. Gore & Associates.

© 2021 W. L. Gore & Associates, Inc. 2198700-EN JUNE 2021

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

