

Remember GORE-TEX® Suture:

The Perfect Close to Your Vascular Procedures

				Thread Size:								
				CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2		
				Approx USP Size:								
				7-0	6-0	5-0	4-0	3-0	2-0	0		
	Double-armed Needles	Needle Code	Ethicon Reference	Syneture Reference	Thread Length	Catalogue Number						
Taper Point 3/8 Circle ⊙		TTc-7	BV130-5	CV-351, MV-135-5	24" – 61 cm	8K10						
		TTc-9	BV-1, BV175-8	CV-1, CVF-1, MV-175-8	18" – 46 cm 24" – 61 cm 30" – 76 cm	8J02 8K02 8M02	7J02 7K02 7M02	6J02 6K02 6M02				
		TTc-12	BV	CV-11, CVF	24" – 61 cm 30" – 76 cm		7M12	6M10				
		TTc-13	C-1	CV-11, CVF-11	18" – 46 cm 24" – 61 cm 30" – 76 cm 36" – 91 cm		7J06 7M06	6J04 6K04 6M04 6N04	5K02 5M02 5N02			
		TT-17	BB	CV-15, CVF-15	36" – 91 cm				5N06	4N04	3N02	
Taper Point 1/2 Circle ⊙		THc-9	N/A	CV-20	30" – 76 cm			6M24				
		THc-13	RB-2	CV-22, CVF-22	18" – 46 cm 24" – 61 cm 30" – 76 cm 36" – 91 cm		7J10 7M10	6J12 6K12 6M12	5K12 5M12 5N12			
		TH-18	RB-1	CV-23, CVF-23	36" – 91 cm 48" – 122 cm			6N16	5N04	4N02	3N04	
		TH-22 (with pledgets)	SH-1	CV-25	36" – 91 cm				3201†† 3202†	3200††		
		TH-26	SH	V-20	36" – 91 cm 48" – 122 cm				5N20	4N06	3N08	2N02
		TH-35	MH	V-26	36" – 91 cm 48" – 122 cm					4N12	3N06	
		TH-35	MH	V-26	36" – 91 cm 48" – 122 cm						3U06	
Piercing Point 3/8 Circle Multifaceted tip for tough tissue, tapered body ⊙		PT-9	CC	KV-1	18" – 46 cm 24" – 61 cm 30" – 76 cm		7J04 8K08 8M08	7K04 7M04	6K06 6M06			
		PT-13	CC-1	KV-11	24" – 61 cm 30" – 76 cm 36" – 91 cm			7K08 7M08	6K08 6M08	5K08 5N08		
		PT-17	V-4	KV-15	36" – 91 cm					5N10	4N10	3N14
Piercing Point 1/2 Circle ⊙		PH-13	N/A	N/A	30" – 76 cm			6M14	5M16			
		PH-17	V-5	KV-5	36" – 91 cm				5N14	4N08	3N12	
		PH-24	V-7	KV-7	36" – 91 cm 48" – 122 cm					4N16	3N10	2N10
		PH-24	V-7	KV-7	36" – 91 cm 48" – 122 cm					4U16	3U10	2U10
		PH-30	N/A	N/A	42" – 107 cm						3P24	
	PH-35	V-34	KV-34	48" – 122 cm							3U20	

Double-armed Needles

† With ePTFE pledgets (3) †† With ePTFE pledgets (4)

Catalogue number in bold type indicates 1:1 needle to thread ratio to minimize suture hole bleeding.

To order a three-dozen box, suffix the catalogue number with the letter B. To order a one-dozen box, suffix the catalogue number with the letter A. Some catalogue numbers available in A (one-dozen) boxes only.

Telephone Orders: 800.528.8763, 928.779.2771 (US) • 00800.6334.4673 (EU)

Needle to Thread Ratio

Thread size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2
Approx USP size:	7-0	6-0	5-0	4-0	3-0	2-0	0
	TTc-7 (1.6:1)	TTc-9 (1.4:1)	TTc-9 (1:1)	TTc-13 (1:1)	TT-17 (1.2:1)	TT-17 (1:1)	TH-26 (1:1)
	TTc-9 (1.9:1)	TTc-12 (1.4:1)	TTc-12 (1:1)	TT-17 (1.2:1)	TH-18 (1.2:1)	TH-18 (1:1)	PH-24 (1:1)
	TTc-13 (2.2:1)	TTc-13 (1.6:1)	TTc-13 (1.4:1)	THc-13 (1:1)	TH-22 (1.3:1)	TH-26 (1.2:1)	
	PT-9 (1.9:1)	THc-13 (1.6:1)	THc-9 (1:1)	TH-18 (1.2:1)	TH-26 (1.6:1)	TH-35 (1.2:1)	
	PT-13 (2.2:1)	PT-9 (1.4:1)	THc-13 (1.4:1)	TH-22 (1.7:1)	TH-35 (1.6:1)	PT-17 (1:1)	
		PT-13 (1.6:1)	TH-18 (1.6:1)	TH-26 (1.7:1)	PT-17 (1.2:1)	PH-17 (1:1)	
			PT-9 (1:1)	PT-13 (1:1)	PH-17 (1.2:1)	PH-24 (1.2:1)	
			PT-13 (1.4:1)	PT-17 (1.2:1)	PH-24 (1.6:1)	PH-30 (1.3:1)	
			PH-13 (1.4:1)	PH-13 (1:1)	PH-35 (1.8:1)	PH-35 (1.5:1)	
				PH-17 (1.2:1)			

Minimize suture hole bleeding with 1:1 needle to thread ratio

Commonly Requested GORE-TEX® Sutures

Proximal Aorto-bifemoral					
Thread Size	Thread Length (cm)	Taper Point Needles	Catalog Number	Piercing Point Needles	Catalog Number
CV-3	91	TH-26	3N08	PH-24	3N10
CV-4	91	TH-26	4N06	PH-24	4N16
	91	TH-18	4N02		
Distal Aorto-bifemoral					
	76	THc-13	5M12		
CV-5	76	TTc-13	5M02		
	91	TH-18	5N04		
CV-6	76	TTc-12	6M10	PT-13	6M08
Chordae Tendinae					
CV-4	91	TH-22	3200*		
	91/122	TH-18	4N02/4U02		
CV-5	91	TH-22	3201*		
	91	TH-18	5N04		
Proximal Femoro-popliteal					
CV-5	91	TTc-13	5N02	PT-13	5N08
CV-6	76	TTc-9	6M02		
CV-6	76	TTc-12	6M10		
CV-6	76	TTc-13	6M04	PT-13	6M08
Distal Femoro-popliteal					
CV-6	76	TTc-9	6M02	PT-9	6M06
CV-7	76	TTc-9	7K02		
	76	TTc-9	7M02	PT-9	7M04
CV-8	76	TTc-9	8K02		

Carotid Endarterectomy					
Thread Size	Thread Length (cm)	Taper Point Needles	Catalog Number	Piercing Point Needles	Catalog Number
CV-6	76	TTc-9	6M02		
	76	TTc-13	6M04		
Axillo-Femoral					
CV-5	91	TTc-13	5N02	PT-13	5N08
	76	TTc-9	6M02		
CV-6	76	TTc-12	6M10		
	76	TTc-13	6M04	PT-13	6M08
AV Access Graft					
CV-6	61/76	TTc-9	6K02/6M02		
	61/76	TTc-12	6K10/6M10		
	61/76	TTc-13	6K04/6M04		
CV-7	61/76	TTc-9	7K02/7M02		
AV Fistulas					
CV-7	46	TTc-9	7J02		
CV-8	46	TTc-9	8J02		

Catalog number in bold type indicates 1:1 needle to thread ratio to minimize suture hole bleeding
* with ePTFE pledgets (4)

For a full listing of available GORE-TEX® Suture, visit goremedical.com.



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)
00800.6334.4673 (Europe) 928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.
GORE®, GORE-TEX®, and designs are trademarks of W. L. Gore & Associates.
Synature is a trademark of Covidien AG or its affiliate.
© 2008, 2009, 2012 W. L. Gore & Associates, Inc. AM3020-EN3 APRIL 2012