


# *Designed to Deliver*


Advanced accessories exclusively  
tailored for TIPS procedures.



PERFORMANCE  
by design

# GORE TIPS Needle

Status band confirms needle tip protection




*Only the GORE TIPS Needle features a status band that when visible confirms that the needle tip is covered.*

Stiff 16 Gauge needle provides responsive needle torque control



Sharp needle tip aids precise needle puncture



Ergonomic handle enhances grip and needle control



Direction indicators verify component alignment during device preparation and initial insertion



TIPS  
NEEDLE

# GORE® TIPS Sheath

Interlocking feature secures dilator hub and 10 Fr introducer sheath during advancement



Direction indicators verify component alignment during device preparation and initial insertion




Transparent tapered hub provides clear view of access sleeve and device delivery catheter




TIPS  
SHEATH



The GORE TIPS Set consists of the GORE® TIPS Sheath and the GORE TIPS Needle, which may be supplied together or individually.



The GORE® TIPS Sheath consists of a flexible 10 Fr introducer sheath with a hemostatic valve and a 10 Fr dilator.



The GORE TIPS Needle consists of a high-strength, stiff 16 Gauge stainless steel hollow shaft needle and 10 Fr guiding catheter.

### ▶ GORE TIPS Set Ordering Information

CONFIGURATIONS	CATALOG NUMBER
GORE TIPS Set	TSET1016
GORE® TIPS Sheath	TSTH1040
10 Fr Introducer Sheath	
10 Fr Dilator	
GORE TIPS Needle	TNDL1656
16 Gauge Needle	
10 Fr Guiding Catheter	



#### W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)  
00800.6334.4673 (Europe)

800.437.8181 (United States)  
928.779.2771 (United States)

[goremedical.com](http://goremedical.com)

Products listed may not be available in all markets.

GORE®, PERFORMANCE BY DESIGN, and designs are trademarks of W. L. Gore & Associates.  
© 2016 W. L. Gore & Associates, Inc. AU0667-EN3 AUGUST 2016