


# ***NOW with Radiopaque Markers for Enhanced Visibility***



## **Summary of changes:**

- Addition of four gold radiopaque markers bonded to the graft at each end of endoprosthesis
- Only low profile (0.018" delivery system), 5–8 mm diameter devices incorporate this change
- Catalogue numbers will change for the 0.018" delivery system devices (PAJ configurations)

## **What has NOT changed:**

- Delivery system and profile unchanged
- 0.035" delivery system devices (PAH configurations) do NOT have radiopaque markers on the endoprosthesis

### **W. L. GORE & ASSOCIATES, INC.**

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)    800.437.8181 (United States)  
00800.6334.4673 (Europe)    928.779.2771 (United States)

[goremedical.com](http://goremedical.com)

The GORE® VIABAHN® Endoprosthesis with PROPATEN Bioactive Surface is sold in some markets as the GORE® VIABAHN® Endoprosthesis with Heparin Bioactive Surface.

Products listed may not be available in all markets.

GORE®, PROPATEN, VIABAHN®, and designs are trademarks of W. L. Gore & Associates.  
© 2017 W. L. Gore & Associates, Inc. AW1772-EN1 NOVEMBER 2017

