




CONFIGURATIONS

GORE® VIATORR® TIPS Endoprosthesis with Controlled Expansion

CATALOGUE NUMBER	ENDOPROSTHESIS INTERNAL DIAMETER (mm)*	GRAFT-LINED LENGTH (cm)	UNLINED LENGTH (cm)	MAXIMUM GUIDEWIRE DIAMETER (in)	INTRODUCER SHEATH SIZE (Fr)	MAXIMUM DILATATION BALLOON DIAMETER (mm)*
PTB8104275	8–10	4	2	≤ 0.035	10	10
PTB8105275	8–10	5	2	≤ 0.035	10	10
PTB8106275	8–10	6	2	≤ 0.035	10	10
PTB8107275	8–10	7	2	≤ 0.035	10	10
PTB8108275	8–10	8	2	≤ 0.035	10	10

* The selected balloon diameter may be 8, 9, or 10 mm, and should not exceed 10 mm. A balloon which can reach an inflation pressure of 10 ATM must be selected, and should be inflated to a minimum of 10 ATM.





CONFIGURATIONS

GORE TIPS Set

CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
TSET1016		

GORE TIPS Set includes introducer sheath, dilator, needle and guiding catheter

GORE® TIPS Sheath

CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
TSTH1040		
10 Fr Introducer Sheath	40	≤ 0.035
10 Fr Dilator	47	≤ 0.035

GORE TIPS Needle

CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
TNDL1656		
16 Gauge Needle	56	≤ 0.035
10 Fr Guiding Catheter	49	≤ 0.035

 Consult Instructions for Use

INTENDED USE: The GORE TIPS Set, GORE® TIPS Sheath and GORE TIPS Needle, are intended to be used together for percutaneous transjugular liver access during diagnostic and interventional procedures in patients undergoing a Transjugular Intrahepatic Portosystemic Shunt (TIPS) procedure.

INDICATIONS FOR USE IN THE US: The GORE® VIATORR® TIPS Endoprosthesis is indicated for use in the *de novo* and revision treatment of portal hypertension and its complications such as variceal bleeding, gastropathy, refractory ascites, and / or hepatic hydrothorax. ^{Rx} Only

Products listed may not be available in all markets.

GORE, VIATORR®, and designs are trademarks of W. L. Gore & Associates.
© 2018 W. L. Gore & Associates, Inc. AX0352-EN1 APRIL 2018



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com