




CONFIGURATIONS

GORE® VIATORR® TIPS Endoprosthesis with Controlled Expansion

CATALOGUE NUMBER	ENDOPROSTHESIS INTERNAL DIAMETER (mm)*	GRAFT-LINED LENGTH (cm)	UNLINED LENGTH (cm)	MAXIMUM GUIDEWIRE DIAMETER (in)	INTRODUCER SHEATH SIZE (Fr)	MAXIMUM DILATATION BALLOON DIAMETER (mm)*
PT8104275	8-10	4	2	≤ 0.035	10	10
PT8105275	8-10	5	2	≤ 0.035	10	10
PT8106275	8-10	6	2	≤ 0.035	10	10
PT8107275	8-10	7	2	≤ 0.035	10	10
PT8108275	8-10	8	2	≤ 0.035	10	10

* The selected balloon diameter may be 8, 9, or 10 mm, and should not exceed 10 mm. A balloon which can reach an inflation pressure of 10 ATM must be selected, and should be inflated to a minimum of 10 ATM.





CONFIGURATIONS

GORE TIPS Set

CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
------------------	-----------------------	---------------------------------

TSET1016

GORE TIPS Set includes introducer sheath, dilator, needle and guiding catheter

GORE® TIPS Sheath

CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
------------------	-----------------------	---------------------------------

TSTH1040

10 Fr Introducer Sheath	40	≤ 0.035
10 Fr Dilator	47	≤ 0.035

GORE TIPS Needle


CATALOGUE NUMBER	EFFECTIVE LENGTH (cm)	ACCEPTS GUIDEWIRE DIAMETER (in)
------------------	-----------------------	---------------------------------

TNDL1656

16 Gauge Needle	56	≤ 0.035
10 Fr Guiding Catheter	49	≤ 0.035

 Consult Instructions for Use

INTENDED USE: The GORE TIPS Set, GORE® TIPS Sheath and GORE TIPS Needle, are intended to be used together for percutaneous transjugular liver access during diagnostic and interventional procedures in patients undergoing a Transjugular Intrahepatic Portosystemic Shunt (TIPS) procedure.

INDICATIONS FOR USE UNDER CE MARK: The GORE® VIATORR® TIPS Endoprosthesis is indicated for use in the treatment of portal hypertension and its complications such as: variceal bleeding refractory to, or intolerant of, conventional therapies, inaccessible varices, gastropathy, refractory ascites, and / or hepatic hydrothorax. Refer to *Instructions for Use* at goremedical.com for a complete description of all contraindications, warnings, precautions and adverse events.  Only

Products listed may not be available in all markets.

GORE, VIATORR®, and designs are trademarks of W. L. Gore & Associates.
© 2018 W. L. Gore & Associates, Inc. AX0586-EN1 APRIL 2018



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com