



GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

CONFIGURATIONS



| Catalogue number | Stent labeled / nominal diameter (mm) | Max. post-dilated stent diameter (mm)* | Crimped stent length (mm) | Catheter length (cm) | Guidewire diameter | Introducer sheath size (Fr) |
|------------------|---------------------------------------|--|---------------------------|----------------------|--------------------|-----------------------------|
| BXA051502B | 5 | 8 | 15 | 135 | 0.035" | 7 |
| BXA051902B | 5 | 8 | 19 | 135 | 0.035" | 7 |
| BXA052902B | 5 | 8 | 29 | 135 | 0.035" | 7 |
| BXA053902B | 5 | 8 | 39 | 135 | 0.035" | 7 |
| BXA055902B | 5 | 8 | 59 | 135 | 0.035" | 7 |
| BXA057902B | 5 | 8 | 79 | 135 | 0.035" | 7 |
| BXA061502B | 6 | 8 | 15 | 135 | 0.035" | 7 |
| BXA061902B | 6 | 8 | 19 | 135 | 0.035" | 7 |
| BXA062902B | 6 | 8 | 29 | 135 | 0.035" | 7 |
| BXA063902B | 6 | 8 | 39 | 135 | 0.035" | 7 |
| BXA065902B | 6 | 8 | 59 | 135 | 0.035" | 7 |
| BXA067902B | 6 | 8 | 79 | 135 | 0.035" | 7 |
| BXA071502B | 7 | 11 | 15 | 135 | 0.035" | 7 |
| BXA071902B | 7 | 11 | 19 | 135 | 0.035" | 7 |
| BXA072902B | 7 | 11 | 29 | 135 | 0.035" | 7 |
| BXA073902B | 7 | 11 | 39 | 135 | 0.035" | 7 |
| BXA075902B | 7 | 11 | 59 | 135 | 0.035" | 7 |
| BXA077902B | 7 | 11 | 79 | 135 | 0.035" | 7 |
| BXA082902B | 8 | 11 | 29 | 135 | 0.035" | 7 |
| BXA083902B | 8 | 11 | 39 | 135 | 0.035" | 7 |
| BXA085902B | 8 | 11 | 59 | 135 | 0.035" | 7 |
| BXA087902B | 8 | 11 | 79 | 135 | 0.035" | 8 |
| BXA092902B | 9 | 13 | 29 | 135 | 0.035" | 8 |
| BXA093902B | 9 | 13 | 39 | 135 | 0.035" | 8 |
| BXA095902B | 9 | 13 | 59 | 135 | 0.035" | 8 |
| BXA097902B | 9 | 13 | 79 | 135 | 0.035" | 8 |
| BXA102902B | 10 | 13 | 29 | 135 | 0.035" | 8 |
| BXA103902B | 10 | 13 | 39 | 135 | 0.035" | 8 |
| BXA105902B | 10 | 13 | 59 | 135 | 0.035" | 8 |
| BXA107902B | 10 | 13 | 79 | 135 | 0.035" | 8 |
| BXA112902B | 11 | 16 | 29 | 135 | 0.035" | 8 |
| BXA113902B | 11 | 16 | 39 | 135 | 0.035" | 8 |
| BXA115902B | 11 | 16 | 59 | 135 | 0.035" | 8 |
| BXA117902B | 11 | 16 | 79 | 135 | 0.035" | 8 |
| BXAL082902B | 8 | 16 | 29 | 135 | 0.035" | 7 |
| BXAL083902B | 8 | 16 | 39 | 135 | 0.035" | 7 |
| BXAL085902B | 8 | 16 | 59 | 135 | 0.035" | 8 |
| BXAL087902B | 8 | 16 | 79 | 135 | 0.035" | 8 |

* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

Products listed may not be available in all markets.

GORE, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.
© 2018–2020 W. L. Gore & Associates, Inc. AX1673-EN2 JANUARY 2020



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific) 800 437 8181 (United States)
00800 6334 4673 (Europe) 928 779 2771 (United States)

goremedical.com