

GORE® CARDIOFORM

ASD Occluder

SIZING

Catalogue number	Treatment range measured with stop flow balloon sizing	Catheter size*
ASD27A	8–15 mm	10 Fr
ASD32A	13–20 mm	10 Fr
ASD37A	18–25 mm	11 Fr
ASD44A	23–30 mm	12 Fr
ASD48A	28–35 mm	14 Fr

* If a 0.035" guidewire is used it is recommended to increase the introducer sheath size by 2 Fr.

Together, improving life



Refer to *Instructions for Use* for a complete description of all warnings, precautions, and contraindications. ^{Rx only}

Products listed may not be available in all markets.

GORE®, CARDIOFORM, TOGETHER, IMPROVING LIFE, and designs are trademarks of W. L. Gore & Associates, Inc.
© 2019 W. L. Gore & Associates, Inc. JUNE 2019

W. L. Gore & Associates, Inc.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific) 800.437.8181 (United States)

00800.6334.4673 (Europe) 928.779.2771 (United States) goremedical.com

