



# GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

## CONFIGURATIONS



| Catalogue number | Stent labeled/<br>nominal<br>diameter (mm) | Max. post-dilated<br>stent diameter<br>(mm)* | Crimped stent<br>length (mm) | Catheter<br>length (cm) | Guidewire<br>diameter | Introducer<br>sheath size<br>(Fr) |
|------------------|--------------------------------------------|----------------------------------------------|------------------------------|-------------------------|-----------------------|-----------------------------------|
| BXA051502W       | 5                                          | 8                                            | 15                           | 135                     | .035"                 | 7                                 |
| BXA051902W       | 5                                          | 8                                            | 19                           | 135                     | .035"                 | 7                                 |
| BXA052902W       | 5                                          | 8                                            | 29                           | 135                     | .035"                 | 7                                 |
| BXA053902W       | 5                                          | 8                                            | 39                           | 135                     | .035"                 | 7                                 |
| BXA055902W       | 5                                          | 8                                            | 59                           | 135                     | .035"                 | 7                                 |
| BXA057902W       | 5                                          | 8                                            | 79                           | 135                     | .035"                 | 7                                 |
| BXA061502W       | 6                                          | 8                                            | 15                           | 135                     | .035"                 | 7                                 |
| BXA061902W       | 6                                          | 8                                            | 19                           | 135                     | .035"                 | 7                                 |
| BXA062902W       | 6                                          | 8                                            | 29                           | 135                     | .035"                 | 7                                 |
| BXA063902W       | 6                                          | 8                                            | 39                           | 135                     | .035"                 | 7                                 |
| BXA065902W       | 6                                          | 8                                            | 59                           | 135                     | .035"                 | 7                                 |
| BXA067902W       | 6                                          | 8                                            | 79                           | 135                     | .035"                 | 7                                 |
| BXA071502W       | 7                                          | 11                                           | 15                           | 135                     | .035"                 | 7                                 |
| BXA071902W       | 7                                          | 11                                           | 19                           | 135                     | .035"                 | 7                                 |
| BXA072902W       | 7                                          | 11                                           | 29                           | 135                     | .035"                 | 7                                 |
| BXA073902W       | 7                                          | 11                                           | 39                           | 135                     | .035"                 | 7                                 |
| BXA075902W       | 7                                          | 11                                           | 59                           | 135                     | .035"                 | 7                                 |
| BXA077902W       | 7                                          | 11                                           | 79                           | 135                     | .035"                 | 7                                 |
| BXA082902W       | 8                                          | 11                                           | 29                           | 135                     | .035"                 | 7                                 |
| BXA083902W       | 8                                          | 11                                           | 39                           | 135                     | .035"                 | 7                                 |
| BXA085902W       | 8                                          | 11                                           | 59                           | 135                     | .035"                 | 7                                 |
| BXA087902W       | 8                                          | 11                                           | 79                           | 135                     | .035"                 | 8                                 |
| BXA092902W       | 9                                          | 13                                           | 29                           | 135                     | .035"                 | 8                                 |
| BXA093902W       | 9                                          | 13                                           | 39                           | 135                     | .035"                 | 8                                 |
| BXA095902W       | 9                                          | 13                                           | 59                           | 135                     | .035"                 | 8                                 |
| BXA097902W       | 9                                          | 13                                           | 79                           | 135                     | .035"                 | 8                                 |
| BXA102902W       | 10                                         | 13                                           | 29                           | 135                     | .035"                 | 8                                 |
| BXA103902W       | 10                                         | 13                                           | 39                           | 135                     | .035"                 | 8                                 |
| BXA105902W       | 10                                         | 13                                           | 59                           | 135                     | .035"                 | 8                                 |
| BXA107902W       | 10                                         | 13                                           | 79                           | 135                     | .035"                 | 8                                 |
| BXA112902W       | 11                                         | 16 <sup>†</sup>                              | 29                           | 135                     | .035"                 | 8                                 |
| BXA113902W       | 11                                         | 16 <sup>†</sup>                              | 39                           | 135                     | .035"                 | 8                                 |
| BXA115902W       | 11                                         | 16 <sup>†</sup>                              | 59                           | 135                     | .035"                 | 8                                 |
| BXA117902W       | 11                                         | 16 <sup>†</sup>                              | 79                           | 135                     | .035"                 | 8                                 |
| BXAL082902W      | 8                                          | 16 <sup>†</sup>                              | 29                           | 135                     | .035"                 | 7                                 |
| BXAL083902W      | 8                                          | 16 <sup>†</sup>                              | 39                           | 135                     | .035"                 | 7                                 |
| BXAL085902W      | 8                                          | 16 <sup>†</sup>                              | 59                           | 135                     | .035"                 | 8                                 |
| BXAL087902W      | 8                                          | 16 <sup>†</sup>                              | 79                           | 135                     | .035"                 | 8                                 |

\* Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

† Technical limit of the device as determined by in-vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the Gore VBX FLEX Clinical Study.

Refer to *Instructions for Use* for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. <sup>Bx only</sup>

Products listed may not be available in all markets.

GORE, VBX, VIABAHN and designs are trademarks of W. L. Gore & Associates.  
© 2021 W. L. Gore & Associates, Inc. 21393952-EN DECEMBER 2021



**W. L. GORE & ASSOCIATES, INC.**  
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific) 0800 6334 4673 (Europe)  
1800 680 424 (Australia/New Zealand) 800 437 8181 (United States) 928 779 2771 (United States)

[goremedical.com](http://goremedical.com)