





PREDICTABLE PERFORMANCE WITH ENHANCEMENTS IN ACCESS AND DELIVERY

Large diameter configurations with lower delivery profile

- Expanded accessibility through smaller sheaths.
- Smaller sheath sizes result in a lower risk of vascular access complications in select patient populations.¹

Endoprosthesis labeled diameter (mm) [†]	Recommended vessel diameter (mm) [‡]	Legacy profile (Fr)	New lower profile (Fr)
9	7.6–8.5	○ 9	○ 8
10	8.6–9.5	○ 11	○ 8
11	9.6–10.5	○ 11	○ 10
13	10.6–12.0	○ 12	○ 10 [§]

- Enhanced visualization under fluoroscopy
 - Four gold radiopaque markers on the distal and proximal ends for more confident delivery in complex cases



* As used by Gore, Heparin Bioactive Surface refers to Gore's proprietary CBAS Heparin Surface.


† Labeled device diameters and lengths are nominal.

‡ Recommended endoprosthesis compression within the vessel is approximately 5–20%.

§ The 13 mm diameter device is not compatible with the 10 Fr COOK FLEXOR CHECK-FLO Introducer.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

Continued innovation for durable outcomes and versatility



Reference

1. Applegate RJ, Sacrinty MT, Kutcher MA, *et al.* Trends in vascular complications after diagnostic cardiac catheterization and percutaneous coronary intervention via the femoral artery, 1998 to 2007. *JACC: Cardiovascular Interventions* 2008;1(3):317–326.

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Only

Products listed may not be available in all markets.

COOK, CHECK-FLO and FLEXOR are trademarks of Cook Medical, Inc.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.

GORE, *Together, improving life*, HEMOBAHN, VIABAHN and designs are trademarks of W. L. Gore & Associates, Inc.

© 2023 W. L. Gore & Associates, Inc. 231290882-EN DECEMBER 2023

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

