

GORE® TAG®

Conformable Thoracic Stent Graft
with ACTIVE CONTROL System

WHERE CONFORMABILITY MEETS CONTROL

Combining the precise control of an intuitive delivery system with the exceptional conformability of the GORE® TAG® Conformable Thoracic Stent Graft, clinicians can achieve predictable patient outcomes.¹⁻³

- Staged deployment enables refinement of placement
- Optional angulation control at intermediate and full diameter
- Hemodynamic stability throughout the deployment
- Expanded large-diameter taper sizes



Together, improving life



GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System

Catalogue number	Intended aortic diameter (mm)	Proximal diameter (mm)	Distal diameter (mm)	Endoprosthesis length (cm)	Recommended GORE® DRYSEAL Flex Introducer Sheath size (Fr)	Oversizing range (%)	Partially uncovered stent length (mm)
Straight Devices							
TGM212110	16-19.5	21	21	10	18	8-31	3
TGM262610	19.5-24	26	26	10	20	8-33	4
TGM282810	22-26	28	28	10	20	8-27	4
TGM282815	22-26	28	28	15	20	8-27	4
TGMR313110	24-29	31	31	10	20	7-29	4
TGMR313115	24-29	31	31	15	20	7-29	4
TGMR313120	24-29	31	31	20	20	7-29	4
TGM343410	27-32	34	34	10	22	6-26	5
TGM343415	27-32	34	34	15	22	6-26	5
TGM343420	27-32	34	34	20	22	6-26	5
TGMR373710	29-34	37	37	10	22	9-28	5
TGMR373715	29-34	37	37	15	22	9-28	5
TGMR373720	29-34	37	37	20	22	9-28	5
TGMR404010	31-37	40	40	10	22	8-29	6
TGMR404015	31-37	40	40	15	22	8-29	6
TGMR404020	31-37	40	40	20	22	8-29	6
TGM454510	34-42	45	45	10	24	7-32	6.5
TGM454515	34-42	45	45	15	24	7-32	6.5
TGM454520	34-42	45	45	20	24	7-32	6.5
Tapered Devices							
TGM262110	19.5-24/16-19.5	26	21	10	20	8-33	4
TGMR312610	24-29/19.5-24	31	26	10	20	7-33	4
TGM342815	27-32/22-26	34	28	15	22	6-27	5
TGMR373115	29-34/24-29	37	31	15	22	7-29	5
TGMR403415	31-37/27-32	40	34	15	22	6-29	6
TGM453715	34-42/29-34	45	37	15	24	7-32	6.5

References

- Farber MA, Krishnasastry KV, Desai N, Starnes BW, Matsumura JS, Tohill BC; TAG 08-02 Clinical Trial Investigators. Five-year outcomes with Conformable GORE TAG Endoprosthesis used in traumatic aortic transections. *Ann Thorac Surg.* 2022;113(5):1536-1542. doi:10.1016/j.athoracsur.2021.05.054
- Jordan WD, Desai N, Letter AJ, Matsumura JS. Long-term outcomes of the Conformable TAG Thoracic Endoprosthesis in a prospective multicenter trial. *J Vasc Surg.* 2021;74(5):1491-1498. doi:10.1016/j.jvs.2021.04.063
- Desai N, Tohill B, Matsumura J, Cambria RP. A multicenter clinical trial on the five-year outcomes following treatment of acute, complicated type B aortic dissection with a conformable stent graft. *Ann Vasc Surg.* 2025;113:1-12. doi:10.1016/j.avsg.2024.12.008

 Consult Instructions for Use at eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. ^{Rx Only}

INDICATIONS FOR USE: The GORE® TAG® Conformable Thoracic Stent Graft is intended for endovascular repair of all lesions of the descending thoracic aorta, including: isolated lesions in patients who have appropriate anatomy, including: adequate iliac/femoral access, aortic inner diameter in the range of 16-42 mm, ≥ 20 mm non-aneurysmal aorta proximal and distal to the lesion; Type B dissections in patients who have appropriate anatomy, including: adequate iliac/femoral access, ≥ 20 mm landing zone proximal to the primary entry tear; proximal extent of the landing zone must not be dissected, diameter at proximal extent of proximal landing zone in the range of 16-42 mm. **CONTRAINDICATIONS:** Patients with known sensitivities or allergies to the device materials; patients who have a condition that threatens to infect the graft.

Products listed may not be available in all markets.

© 2026 W. L. Gore & Associates, Inc. All rights reserved. All trademarks referenced are trademarks of either a member of the Gore group of affiliated companies or their respective owners. "Together, improving life" mark and designs are trademarks of a Gore company.
25AR1004-EN02 MARCH 2026 MAT-0244-1

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

