




# GORE® VIATORR® TIPS 覆膜支架系

器械规格

# GORE® VIATORR® TIPS 覆膜支架系统

专为门静脉高压及其并发症的初治和修复治疗设计

| 产品编号       | 标示内径 (mm) | 覆膜段长度 (cm) | 裸段长度 (cm) | 最大导丝直径 (英寸) | 止血导引鞘 (Fr) | 最大扩张球囊直径 (mm) |
|------------|-----------|------------|-----------|-------------|------------|---------------|
| PTB084275W | 8         | 4          | 2         | ≤ 0.038     | 10         | 8             |
| PTB085275W | 8         | 5          | 2         | ≤ 0.038     | 10         | 8             |
| PTB086275W | 8         | 6          | 2         | ≤ 0.038     | 10         | 8             |
| PTB087275W | 8         | 7          | 2         | ≤ 0.038     | 10         | 8             |
| PTB088275W | 8         | 8          | 2         | ≤ 0.038     | 10         | 8             |
| PTB104275W | 10        | 4          | 2         | ≤ 0.038     | 10         | 10            |
| PTB105275W | 10        | 5          | 2         | ≤ 0.038     | 10         | 10            |
| PTB106275W | 10        | 6          | 2         | ≤ 0.038     | 10         | 10            |
| PTB107275W | 10        | 7          | 2         | ≤ 0.038     | 10         | 10            |
| PTB108275W | 10        | 8          | 2         | ≤ 0.038     | 10         | 10            |

**适用范围:** GORE® VIATORR® TIPS 覆膜支架适用于初治和修复治疗门静脉高血压及其并发症,如静脉曲张出血、胃病、顽固性腹水和 / 或肝性胸水。

 Consult Instructions  
for Use  
[eifu.goremedical.com](http://eifu.goremedical.com)

请参阅 [eifu.goremedical.com](http://eifu.goremedical.com) 上的“使用说明”,以获取该产品销售市场的所有适用的适应症、警告、注意事项和禁忌症的完整说明。IX Only

所列产品可能并非在所有市场中都有销售。

GORE, VIATORR and designs are trademarks of W. L. Gore & Associates.  
© 2020 W. L. Gore & Associates, Inc. AZ0872-SZH1 APRIL 2020



**W. L. GORE & ASSOCIATES, INC.**  
Flagstaff, AZ 86004

+65 67332882 (Asia Pacific)  
1800 680 424 (Australia / New Zealand)  
00800 6334 4673 (Europe)  
800 437 8181 (United States)  
928 779 2771 (United States)

[goremedical.com](http://goremedical.com)