

GORE® EXCLUDER® Thoracoabdominal Branch Endoprosthesis (TAMBE) Selection Form

Patient confidential information — Do not disclose legally protected data.

The following information is required to ensure that the appropriate devices and backups are available for the procedure.

Patient ID:

Physician:

Institution:

Date:

Aortic Component: Proximal diameter (mm) x distal diameter (mm) x length (mm)					
Qty.	Catalogue number	Device size	Intended proximal aortic inner diameter (mm)	Recommended sheath size (Fr)	
	ATAA43120160A	31 x 20 x 160	22 – 29	22	
	ATAA43720160A	37 x 20 x 160	27 – 34	22	
Distal Bifurcated Component: Proximal diameter (mm) x distal diameter (mm) x length (cm)					
Qty.	Catalogue number	Device size	Intended iliac vessel diameter (mm)	Recommended sheath size (Fr)	
	CEB231010A	23 x 10 x 10	8.0 – 9.0	16	
	CEB231210A	23 x 12 x 10	10.0 – 11.0	16	
	CEB231410A	23 x 14.5 x 10	12.0 – 13.5	16	
Branch components (GORE® VIABAHN® VBx Balloon Expandable Endoprosthesis): Diameter (mm) x length (mm)					
Initial device					
Qty.	Catalogue number	Nominal device size	Final stent ID range	Catalogue number	Qty.
	BXA053902A	5 x 39	4.5 – 8.0	BXB053902A	
	BXA055902A	5 x 59	4.5 – 8.0	BXB055902A	
	BXA057902A	5 x 79	4.3 – 8.0	BXB057902A	
	BXA063902A	6 x 39	5.0 – 8.0	BXB063902A	
	BXA065902A	6 x 59	5.1 – 8.0	BXB065902A	
	BXA067902A	6 x 79	5.3 – 8.0	BXB067902A	
	BXA073902A	7 x 39	6.0 – 11.0	BXB073902A	
	BXA075902A	7 x 59	6.0 – 11.0	BXB075902A	
	BXA077902A	7 x 79	6.0 – 11.0	BXB077902A	
	BXA083902A	8 x 39	6.7 – 11.0	BXB083902A	
	BXA085902A	8 x 59	6.8 – 11.0	BXB085902A	
	BXA087902A	8 x 79	7.0 – 11.0	BXB087902A	
	BXA093902A	9 x 39	8.0 – 13.0	BXB093902A	
	BXA095902A	9 x 59	7.8 – 13.0	BXB095902A	
	BXA097902A	9 x 79	7.8 – 13.0	BXB097902A	

Accessories

GORE® DRYSEAL Flex Introducer Sheath			
Qty.	Catalogue number*	Sheath size (Fr)	Nominal sheath OD (mm)
	GDSFI233	12	4.7
	GDSFI245	12	4.7
	GDSFI433	14	5.3
	GDSFI533	15	5.6
	GDSF2233	22	8.2
GORE® Tri-Lumen Catheter			
Qty.	Catalogue number	Device length	Size (Fr)
	TLC140	140 cm	7.5
GORE® Tri-Lobe Balloon Catheter			
Qty.	Catalogue number	Size (Fr)	Vessel diameter (mm)
	BCM1634	18	16 – 32
	BCL2645	18	26 – 42
GORE® Molding & Occlusion Balloon			
Qty.	Catalogue number	Size (Fr)	Balloon diameter (mm)
	MOB37	10	10 – 37

* GDSF catalogue numbers are not available in all regions; use DSF catalogue numbers instead where appropriate.

Contralateral Leg Endoprosthesis (GORE® EXCLUDER® AAA Endoprosthesis): Distal diameter (mm) x length (cm)				
Qty.	Catalogue number	Device size*	Intended iliac vessel diameter (mm)	Recommended sheath size (Fr)
	PLC121000	12 x 10	10 – 11	12
	PLC121200	12 x 12	10 – 11	12
	PLC121400	12 x 14	10 – 11	12
	PLC141000	14.5 x 10	12 – 13.5	12
	PLC141200	14.5 x 12	12 – 13.5	12
	PLC141400	14.5 x 14	12 – 13.5	12
	PLC161000*	16 x 9.5	13.5 – 14.5	12
	PLC161200*	16 x 11.5	13.5 – 14.5	12
	PLC161400*	16 x 13.5	13.5 – 14.5	12
	PLC181000*	18 x 9.5	14.5 – 16.5	12
	PLC181200*	18 x 11.5	14.5 – 16.5	12
	PLC181400*	18 x 13.5	14.5 – 16.5	12
	PLC201000*	20 x 9.5	16.5 – 18.5	12
	PLC201200*	20 x 11.5	16.5 – 18.5	12
	PLC201400*	20 x 13.5	16.5 – 18.5	12
	PLC231000*	23 x 10	18.5 – 21.5	14
	PLC231200*	23 x 12	18.5 – 21.5	14
	PLC231400*	23 x 14	18.5 – 21.5	14
	PLC271000*	27 x 10	21.5 – 25.0	15
	PLC271200*	27 x 12	21.5 – 25.0	15
	PLC271400*	27 x 14	21.5 – 25.0	15
*PLEASE NOTE: All large diameter Contralateral Leg Endoprostheses (16, 18, 20, 23, 27 mm) can be used as Iliac Extenders.				
Iliac Extender (GORE® EXCLUDER® Iliac Branch Endoprosthesis): Distal diameter (mm) x length (cm)				
Qty.	Catalogue number	Device size	Intended iliac vessel diameter (mm)	Recommended sheath size (Fr)
	PLL161007	10 x 7	8.0 – 9.0	12
	PLL161207	12 x 7	10.0 – 11.0	12
	PLL161407	14.5 x 7	12.0 – 13.5	12

Additional devices

Distal Bifurcated Component Extender (GORE® EXCLUDER® AAA Endoprosthesis) Distal diameter (mm) x length (cm)				
Qty.	Catalogue number	Device size	Intended iliac vessel diameter (mm)	Recommended sheath size (Fr)
	PLA230300	23 x 3.3	19 – 21	16
Other items				
Qty.	Catalogue number	Description		

Products listed may not be available in all markets.

GORE, ACTIVE CONTROL, DRYSEAL, EXCLUDER, VIABAHN and designs are trademarks of W. L. Gore & Associates. © 2023, 2025 W. L. Gore & Associates, Inc. 24ARI240-EN02 MAY 2025

