

GORE® VIABAHN® VBX

Balloon Expandable Endoprosthesis

6 Fr
compatible

GREATER
VERSATILITY



COMPETITIVE CONFIGURATIONS COMPARISON CARD

Together, improving life



The VBX Stent Graft offers the **most** 6 Fr compatible configurations and the **longest** device length among balloon expandable stent grafts.¹⁻⁴

	GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis ^a length (mm)						ATRIUM® iCAST® Covered Stent length (mm) ^b					BD® LIFESTREAM™ Balloon Expandable Vascular Covered Stent length (mm)			
	15	19	29	39	59	79	16	22	29 ^c /32	38/41 ^c	59/61 ^c	16	26	37/38	58
5	6 Fr						6 Fr		7 Fr						
6															
7															
8	7 Fr											7 Fr			
8L															
9															
10	8 Fr											8 Fr			
11															
12															
									9 Fr						

a. Latest VBX Stent Graft profile (Fr).

b. If using a 10 mm device with a COOK® Introducer Sheath, it is recommended to use an 8 Fr sheath to minimize the potential for resistance.

c. 29, 41, and 61 mm lengths on iCAST® Covered Stent only apply to the 12 mm diameter configuration.

The VBX Stent Graft offers the **widest** range of diameter adjustability among balloon expandable stent grafts.¹⁻⁴

Nominal diameter (mm)	Max post-dilated stent diameter (mm)		
	GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis ^a	ATRIUM® iCAST® Covered Stent	BD® LIFESTREAM™ Balloon Expandable Vascular Covered Stent
5	8	7.3-9.8	10
6	8	7.3-10	10
7	11	7.3-10.1	10
8	11	9.3-10.2	10
8L	16 ^b	NA	NA
9	13	9.3-10.4	12
10	13	10.6	12
11	16 ^b	NA	NA
12	NA	14.9-15.1 ^c	12

Abbreviation: NA, not applicable.

All VBX Stent Graft configurations available in 0.035" over-the-wire 80 and 135 cm catheter lengths.

a. Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

b. Technical limit of the device as determined by in vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the Gore VBX FLEX Clinical Study. See *Instructions for Use*.

c. Device expansion beyond 12 mm is outside the reference vessel diameters listed in the indication statement of the iCAST® Covered Stent.

1. W. L. Gore & Associates. *GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis Instructions for Use*. W. L. Gore & Associates; 2025. MD169334, Rev 8. Accessed April 28, 2026. <https://eifu.goremedical.com/>
2. BD. *The Control You Need to Deliver Accurate Treatment: LIFESTREAM™ Balloon Expandable Covered Stent*. BD; 2022. BD-52176v2. Accessed June 9, 2026. <https://www.bd.com/content/dam/bd-assets/na/peripheral-intervention/web-assets/us/documents/bd-52176-lifestream-br.pdf>
3. Atrium Medical Corporation. *iCAST Covered Stent System Instructions for Use*. Atrium Medical Corporation; 2025. EL012022, Rev AA.
4. U.S. Food and Drug Administration. *Summary of Safety and Effectiveness Data (SSED): LifeStream™ Balloon Expandable Vascular Covered Stent*. PMA P160024. Accessed June 2, 2026. https://www.accessdata.fda.gov/cdrh_docs/pdf16/P160024B.pdf



Consult Instructions for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. **Rx Only**

INDICATIONS FOR USE IN THE US: The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis is indicated for the treatment of de novo or restenotic lesions found in iliac arteries with reference vessel diameters ranging from 5 mm–13 mm and lesion lengths up to 110 mm, including lesions at the aortic bifurcation. The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis is also indicated for use with thoracoabdominal and pararenal branched devices indicated with the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis as a branch component. **CONTRAINDICATIONS:** Do not use the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis in patients with known hypersensitivity to heparin, including those patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II.

Products listed may not be available in all markets.

© 2026 W. L. Gore & Associates, Inc. All rights reserved. All trademarks referenced are trademarks of either a member of the Gore group of affiliated companies or their respective owners. “Together, improving life” mark and designs are trademarks of a Gore company. 26PL1095-EN02 AUGUST 2026 MAT-0596-2

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 **Australia/New Zealand** 1800 680 424
Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

