


GORE® VIABAHN® VBX
Balloon Expandable Endoprosthesis



6Fr
compatible

**GREATER
VERSATILITY**

CONFIGURATIONS

Together, improving life



Catalogue numbers

GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis (VBX Stent Graft)

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm) ^a	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Reduced profile VBX Stent Graft introducer sheath size (Fr)
BXA051501A	BXB051501A	5	8	15	80	.035	6
BXA051901A	BXB051901A	5	8	19	80	.035	6
BXA052901A	BXB052901A	5	8	29	80	.035	6
BXA053901A	BXB053901A	5	8	39	80	.035	6
BXA055901A	BXB055901A	5	8	59	80	.035	6
BXA057901A	BXB057901A	5	8	79	80	.035	6
BXA061501A	BXB061501A	6	8	15	80	.035	6
BXA061901A	BXB061901A	6	8	19	80	.035	6
BXA062901A	BXB062901A	6	8	29	80	.035	6
BXA063901A	BXB063901A	6	8	39	80	.035	6
BXA065901A	BXB065901A	6	8	59	80	.035	6
BXA067901A	BXB067901A	6	8	79	80	.035	6
BXA071501A	BXB071501A	7	11	15	80	.035	6
BXA071901A	BXB071901A	7	11	19	80	.035	6
BXA072901A	BXB072901A	7	11	29	80	.035	6
BXA073901A	BXB073901A	7	11	39	80	.035	6
BXA075901A	BXB075901A	7	11	59	80	.035	6
BXA077901A	BXB077901A	7	11	79	80	.035	6
BXA082901A	BXB082901A	8	11	29	80	.035	7
BXA083901A	BXB083901A	8	11	39	80	.035	7
BXA085901A	BXB085901A	8	11	59	80	.035	7
BXA087901A	BXB087901A	8	11	79	80	.035	7

^a Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm) ^a	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Reduced profile VBX Stent Graft introducer sheath size (Fr)
BXA092901A	BXB092901A	9	13	29	80	.035	7
BXA093901A	BXB093901A	9	13	39	80	.035	7
BXA095901A	BXB095901A	9	13	59	80	.035	7
BXA097901A	BXB097901A	9	13	79	80	.035	7
BXA102901A	BXB102901A	10	13	29	80	.035	7
BXA103901A	BXB103901A	10	13	39	80	.035	7
BXA105901A	BXB105901A	10	13	59	80	.035	7
BXA107901A	BXB107901A	10	13	79	80	.035	7
BXA112901A	BXB112901A	11	16 ^b	29	80	.035	8
BXA113901A	BXB113901A	11	16 ^b	39	80	.035	8
BXA115901A	BXB115901A	11	16 ^b	59	80	.035	8
BXA117901A	BXB117901A	11	16 ^b	79	80	.035	8
BXAL082901A	BXBL082901A	8	16 ^b	29	80	.035	7
BXAL083901A	BXBL083901A	8	16 ^b	39	80	.035	7
BXAL085901A	BXBL085901A	8	16 ^b	59	80	.035	8
BXAL087901A	BXBL087901A	8	16 ^b	79	80	.035	8
BXA051502A	BXB051502A	5	8	15	135	.035	6
BXA051902A	BXB051902A	5	8	19	135	.035	6
BXA052902A	BXB052902A	5	8	29	135	.035	6
BXA053902A	BXB053902A	5	8	39	135	.035	6
BXA055902A	BXB055902A	5	8	59	135	.035	6
BXA057902A	BXB057902A	5	8	79	135	.035	6
BXA061502A	BXB061502A	6	8	15	135	.035	6
BXA061902A	BXB061902A	6	8	19	135	.035	6
BXA062902A	BXB062902A	6	8	29	135	.035	6
BXA063902A	BXB063902A	6	8	39	135	.035	6
BXA065902A	BXB065902A	6	8	59	135	.035	6
BXA067902A	BXB067902A	6	8	79	135	.035	6

^a Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

^b Technical limit of the device as determined by in vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the Gore VBX FLEX Clinical Study. See *Instructions for Use*.

VBX Stent Graft catalogue number	Reduced profile VBX Stent Graft catalogue number	Stent labeled/nominal diameter (mm)	Maximum post-dilated stent diameter (mm) ^a	Stent length (mm)	Catheter length (cm)	Guidewire diameter (in)	Reduced profile VBX Stent Graft introducer sheath size (Fr)
BXA071502A	BXB071502A	7	11	15	135	.035	6
BXA071902A	BXB071902A	7	11	19	135	.035	6
BXA072902A	BXB072902A	7	11	29	135	.035	6
BXA073902A	BXB073902A	7	11	39	135	.035	6
BXA075902A	BXB075902A	7	11	59	135	.035	6
BXA077902A	BXB077902A	7	11	79	135	.035	6
BXA082902A	BXB082902A	8	11	29	135	.035	7
BXA083902A	BXB083902A	8	11	39	135	.035	7
BXA085902A	BXB085902A	8	11	59	135	.035	7
BXA087902A	BXB087902A	8	11	79	135	.035	7
BXA092902A	BXB092902A	9	13	29	135	.035	7
BXA093902A	BXB093902A	9	13	39	135	.035	7
BXA095902A	BXB095902A	9	13	59	135	.035	7
BXA097902A	BXB097902A	9	13	79	135	.035	7
BXA102902A	BXB102902A	10	13	29	135	.035	7
BXA103902A	BXB103902A	10	13	39	135	.035	7
BXA105902A	BXB105902A	10	13	59	135	.035	7
BXA107902A	BXB107902A	10	13	79	135	.035	7
BXA112902A	BXB112902A	11	16 ^b	29	135	.035	8
BXA113902A	BXB113902A	11	16 ^b	39	135	.035	8
BXA115902A	BXB115902A	11	16 ^b	59	135	.035	8
BXA117902A	BXB117902A	11	16 ^b	79	135	.035	8
BXAL082902A	BXBL082902A	8	16 ^b	29	135	.035	7
BXAL083902A	BXBL083902A	8	16 ^b	39	135	.035	7
BXAL085902A	BXBL085902A	8	16 ^b	59	135	.035	8
BXAL087902A	BXBL087902A	8	16 ^b	79	135	.035	8

^a Secondary balloon required to post-dilate the stent beyond its nominal deployed diameter (secondary balloon not included).

^b Technical limit of the device as determined by in vitro testing for the indicated use; device expansion beyond 13 mm was not studied as part of the Gore VBX FLEX Clinical Study. See *Instructions for Use*.



Consult Instructions for Use
eifu.goremedical.com

Refer to *Instructions for Use* at eifu.goremedical.com for a complete description of all applicable indications, warnings, precautions and contraindications for the markets where this product is available. Rx Only

INDICATIONS FOR USE IN THE US: The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis is indicated for the treatment of de novo or restenotic lesions found in iliac arteries with reference vessel diameters ranging from 5 mm–13 mm and lesion lengths up to 110 mm, including lesions at the aortic bifurcation. **CONTRAINDICATIONS:** Do not use the GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis in patients with known hypersensitivity to heparin, including those patients who have had a previous incident of Heparin-Induced Thrombocytopenia (HIT) type II. Products listed may not be available in all markets.

© 2026 W. L. Gore & Associates, Inc. All rights reserved. All trademarks referenced are trademarks of either a member of the Gore group of affiliated companies or their respective owners. “Together, improving life” mark and designs are trademarks of a Gore company. 26PL6022-EN01 AUGUST 2026 MAT-1184-1

W. L. Gore & Associates, Inc.
goremedical.com

Asia Pacific +65 67332882 Australia/New Zealand 1800 680 424 Europe 00800 6334 4673
United States Flagstaff, AZ 86004 800 437 8181 928 779 2771

